# The heart of precision


Nobel Prize®-winning principle behind TRACE technology guarantees highly sensitive and specific test results

Homogeneous assay design ensures excellent analytical precision and high reproducibility by eliminating the need of washing and separation steps.

TRACE™ technology provides simultaneous measurement of signals of bound and unbound antibodies, which reduces effects of possible interferences. The unspecific signals are excluded by time resolved measurement.

## TRACE: Time Resolved Amplified Cryptate Emission

The unique measurement principle underlying the precision and usability of the B-R-A-H-M-S KRYPTOR analyzers

**Amplified cryptate emission:** Assay antibodies are labeled with cryptate, acting as a donor, and fluorophore, acting as energy acceptor. When the target antigen is present in the sample, antibodies build an immunocomplex around the target. Bringing cryptate close to fluorophore results in energy transfer and emission of the signal with the wave length much longer than the ones of unbound antibodies.

Watch an animated visual presentation of TRACE technology  
brahms-instruments.com/technology

## Ease of use and confidence in results

Accurate and precise


Fast time to result

Less waste and reagent consumption

Biotin independent assay design

**Time resolved:** Short signals of unbound fluorophores are differentiated by a 50  $\mu$ s delay of reading window from the long signal of the bound immunocomplex. This long-lived signal is proportional to the concentration of the analyte of interest.

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