

Labclass No.2 – Cleavage of plasmid DNA using restrictionendonucleases

Workflow:

- 3 party: **A) Restriction cleavage of plasmid DNA**
B) Gel preparation (0.8% agarose gel)
C) Plasmid DNA electrophoresis

A) Restriction digestion of plasmid DNA

Method:

- 1) Incubate for **30 minutes**
- 2) **thermal inactivation** of the enzyme - 65 ° C, 5 minutes
- 3) the so-called **L-form** of the plasmid is formed (linearized)

Reaction mixture:

- plasmid DNA **15 µL** (recalculate the actual amount of DNA according to the determined concentration)
 - NEB4 10X **2 µL** buffer
 - BSA 10X **2 µL**
- => centrifugation (short spin, 6000g)
- **Xmn I** **1 µL**
- TOTAL: **20 µL**

B) Gel preparation (0.8% agarose gel)

- we prepare 0.8% agarose gel in TBE buffer.

Method:

- 1) - first it is necessary to glue the pads (into which we pour the gel) and insert combs into them
 - 2) - weigh agarose and pour 100 ml of 0.5% TBE buffer
 - 3) - boil agarose in the microwave (note: the solution must be transparent, 3-4 times boiling, no bubbles)
 - 4) - then the boiled agarose is cooled to 50 ° C under running cold water
 - 5) - add MIDORI GREEN 1:10 000 (4 µl per 100 ml gel)
 - 6) - pour the gel into the pad and let it solidify for 30 minutes.
- Pour 0.5% TBE buffer into the electrophoresis baths.

C) Plasmid DNA electrophoresis

- is performed to control restriction cleavage
- agarose, after solidification, forms a gel = molecular sieve, in which DNA moves (has a negative charge), i.e. it moves from a negatively charged pole (electrode)

Method:

- 1) use 0.8% gel prepared in block B
- 2) remove the comb => this will create holes, remove the adhesive TAPE !!
- 3) insert the gel into the tub so that the wells are above the red stripe
- 4) pour buffer to the MAX line

5) sample pipetting: ladder 3 μ l (1kbp); plasmid DNA: 7 μ l and 2 μ l bromophenol blue - make drops on parafilm, mix

6) each student applies his undigested and digested plasmid to the gel, first always a sample of undigested DNA, then digested; (students - write down the hole numbers with your sample !!!)

7) 150 V, 30 min

After electrophoresis, the gel is transferred to a transilluminator and photographed (under UV light).