Name: Topic:

Micropropagation of Saintpaulia ionantha Wendl.

Transfer of the culture ("subculture")

Plant material: organogenic culture of Saintpaulia ionantha Wendl.

Medium: M-S, sucrose 20 g.1⁻¹, 0.65% agar, pH 5.5

Procedure:

- 1. Transfer tissue culture vessels containing African violet cultures from the culture room into a laminar flow hood.
- 2. Prepare the vessels containing fresh medium and sterile tools (forceps, scalpels), sterile Petri dish.
- 3. Pass the opening of the vessel through the flame.
- 4. Transfer the explants to a sterile Petri dish.
- 5. Remove the necrotic basal tissue of the explants with lower leaves and gently push the isolated apical parts of explants into the fresh medium.
- 6. Pass the opening of the vessel through the flame carefully. Cap the vessels and label.
- 7. Record the date and the number of the culture.
- 8. Culture in the culture room under the light (cool white fluorescent light, photoperiod 16/8 hours, PAR 30 μmol.m⁻².sec⁻¹) at 22°C for 4 weeks.

Evaluation:

Control the contamination of the culture in the following cultivation period.

Notes:

If the culture will be transferred to non-sterile conditions, root the shoots on rooting medium containing lower concentrations of salts (1/3 MS).