

Atomic Layer Deposition (ALD)

Instrument description

Atomic Layer Deposition is a deposition technique for very thin layers with the thickness control down to a single atomic layer. It belongs to the CVD techniques family. The thickness precision is achieved by pulsed deposition, where first a metal-containing precursor is introduced into the chamber and after a short time (allowing for a monolayer adsorption) the chamber is pumped down. Following step is an exposure to the oxidizing precursor (for oxides) or nitrogen containing precursor (for nitrides). Thus, a monolayer of target material is grown. The metal-containing precursors are usually organometallic ones, for oxidation a water or oxygen plasma can be used, nitridation is done using water or nitrogen plasma. To achieve the deposition in the ALD mode, sample is heated up to a certain temperature, for most processes being in the range 150 °C - 300 °C.



Instrument: Ultratech-Cambridge Nanotech Fiji 200

ALD system for up to 8" samples, equipped with plasma generator.
 Standard materials: Al_2O_3 , AlN, HfO_2 , HfN, TiO_2 , TiN, SiO_2 , SiN, other materials on request.


Features:

- thermal deposition within range RT-500 °C
- 4 precursor lines, with possible upgrade to 6
- plasma-enhanced deposition (3 plasma gas lines)
- expo mode for homogeneous deposition on high-aspect-ratio nanostructures
- controlling software allows preparation/ /modification/storage of individual recipes
- fully automatic programmable operation


Application

Al_2O_3 diffusion barrier deposition

XPS analysis



AFM analysis



Technical specification

General info:

Up to 8 inch sample size
Deposition uniformity 1.5 % (1 σ) for Al₂O₃

Delivery/Carrier system:

- solid, liquid and gas precursors heated up to 200 °C
- ALD Booster™ for low pressure precursors
- 3 x mass flow controlled plasma lines, 200 sccm

purge/plasma gas lines

4 precursors lines

control computer

electronics

Operational modes:

- Continuous™ mode
- Exposure™ mode
- Plasma™ mode

plasma chamber/generator

reactor chamber

with sample plate heated up to 500 °C

Deposition modes:

- high speed
- high aspect ratio
- plasma assisted

precursor trap

stop valve

dry pumping system

sample process:

Al₂O₃ from water and TMA

● aluminium
● carbon
● hydrogen
● oxygen

initial surface

ALD CYCLES

precursor flows in

precursor adsorption and desorption of reaction products

chamber purge

second precursor flows in

reaction with the adsorbed layer and desorption of reaction products

chamber purge

ALD CYCLES

thin layer after deposition

Contact

Core Facility: Nanofabrication and Nanocharacterization

Section: Etching & Deposition

Contact person: David Škoda

core.facility@ceitec.vutbr.cz

Detailed information: www.ceitec.eu

Instrument location: Brno University of Technology
Faculty of Mechanical Engineering
Department of Physical Engineering
Technická 2896/2, 616 69 Brno
Czech Republic