

Klonování genu pro glukosa oxidasu z *A. niger*

>AF234246.2:486-2303 *Aspergillus niger* glucose oxidase (GO2) gene, complete cds

ATGCAGACTCTCCTTGTGAGCTCGCTTGTGGTCTCCCTCGCTGCGGCCCTGCCACACTACATCAGGAGCA
ATGGCATTGAAGCCAGCCTCCTGACTGATCCCAAGGATGTCTCCGGCCGCACGGTCGACTACATCATCGC
TGGTGGAGGTCTGACTGGACTCACCACCGCTGCTCGTCTGACGGAGAACCCCAACATCAGTGTGCTCGTC
ATCGAAAGTGGCTCCTACGAGTCGGACAGAGGTCTATCATTGAGGACCTGAACGCCTACGGCGACATCT
TTGGCAGCAGTGTAGACCACGCCTACGAGACCGTGGAGCTCGCTACCAACAATCAAACCGCGCTGATCCG
CTCCGAAATGGTCTCGGTGGCTCTACTCTAGTGAATGGTGGCACCTGGACTCGCCCCACAAGGCACAG
GTTGACTCTTGGGAGACTGTCTTTGAAATGAGGGCTGGAAGTGGGACAATGTGGCCGCTACTCCCTCC
AGGCTGAGCGTGTCTCGCGACCAAATGCCAAACAGATCGCTGCTGGCCACTACTTCAACGCATCCTGCCA
TGGTGTAAATGGTACTGTCCATGCCGGACCCCGCGACACCGGCGATGACTATTCTCCCATCGTCAAGGCT
CTCATGAGCGTGTCTGAAGACCGGGGCGTTCCCAACAAAGAAAGACTTCGGATGCGGTGACCCCATGGTG
TGTCCATGTTCCCAACACCTTGCACGAAGACCAAGTGCCTCCGATGCCGCTCGCGAATGGCTACTTCC
CAACTACCAACGTCCCAACCTGCAAGTCTGACCGGACAGTATGTTGGTAAGGTGCTCCTTAGCCAGAAC
GGCACCACCCCTCGTGCCGTTGGCGTGAATTCGGCACCCACAAGGGCAACACCCACAACGTTTACGCTA
AGCACGAGGTCTCCTGGCCGCGGGCTCCGCTGTCTCTCCACAATCCTCGAATATTCCGGTATCGGAAT
GAAGTCCATCCTGGAGCCCTTGGTATCGACACCGTCTGACCTGCCCGTCCGGCTTGAACCTGCAGGAC
CAGACCACCGCTACCGTCCGCTCCCGCATCACCTCTGCTGGTGCAGGACAGGGACAGGCCGCTTGGTTCG
CCACCTTCAACGAGACCTTTGGTGACTATTCCGAAAAGGCACACGAGCTGCTCAACACCAAGCTGGAGCA
GTGGGCCGAAGAGGCCGCTCGCCCGTGGCGGATTCCACAACACCACCGCCTTGCTCATCCAGTACGAGAAC
TACCGCGACTGGATTGTCAATCACAACGTGCGCTACTCGGAACTCTTCCCTCGACACTGCCGGAGTGGCCA
GCTTCGATGTGTGGGACCTTCTGCCCTTCGACCGAGGATACGTCCACATCCTCGACAAGGACCCCTACCT
CCACCACCTTTCCTACGACCCTCAGTACTTCTCAACGAGCTCGACCTGCTCGGTGAGGCTGCCGCTACT
CAGCTGGCCCGCAACATCTCCAACCTCCGGTGCCATGCAGACCTACTTCCGCTGGGGAGATACTCCCCGGTG
ATAACCTCGCGTATGATGCCGATTTGAGCGCCTGGACTGAGTACATCCCGTACCCTTCCGTCTTAACCTA
CCATGACGTAGGTACTTGTCCATGATGCCGAAGGAGATGGGCAGTGTGTTGATAATGCTGCCCGTGTG
TATGGTGTGCGGGGACTGCGTGTGATTGATGGTCTATTCTCCTACGCAAATGTCGTCCCATGTCATGA
CGGTGTTCTATGCCATGGCTTTGAAAATTTCCGGATGCTATCTTGAAGATTATGCTTCCATGCAG**TGA**

N-HIS,T7-GO2-HIS

*Bam*H1 (GGATCC)

*Xho*I (CTCGAG)

5` TCGCGGATCCCTGCCACACTACATCAGGAGCA

5` GGTGCTCGAGCTGCATGGAAGCATAATCTTCCA

N-HIS,T7-GO2

*Bam*H1 (GGATCC)

*Xho*I (CTCGAG)

5` TCGCGGATCCCTGCCACACTACATCAGGAGCA

5` GGTGCTCGAGTCACTGCATGGAAGCATAATCTTCC

N-HIS-GO2

*Nde*I (CATATG)



*Xho*I (CTCGAG)

5` TATGCATATGCTGCCACACTACATCAGGAGCA

5` GGTGCTCGAGTCACTGCATGGAAGCATAATCTTCC

>AAF59929.2 glucose oxidase [Aspergillus niger]

MQTLLVSSLVSLAAALPHYIRSNGLIEASLLTDPKDVSGRTVDYIIAGGGTGLTTAARLTENPNISVLV
 IESGSYESDRGPIIEDLNAYGDI FGSSVDHAYETVELATNNQTALIRSGNGLGGSTLVNNGGTWTRPHKAQ
 VDSWETVFGNEGWNWDNVAAYS LQAERARAPNAKQIAAGHYFNASCHGVNGTVHAGPRDTGDDYSPIVKA
 LMSAVEDRGVPTKKDFGCGDPHGVSMPFNTLHEDQVRS DAAREWLLPNYQRPNLQVLTGQYVGKVL L SQN
 GTTPRAVGVEFGTHKGNTHNVYAKHEVLLAAGSAVSPTILEYSGIGMKSILEPLGIDTVVDLPVGLNLQD
 QTTATVRSRITSAGAGQGQAAWFATFNETFGDYSEKAHELLNTKLEQWAEAVARGGFHNTTALLIQYEN
 YRDWIVNHNVAyselFLDTAGVASFDVWDL L PDRGYVHILDKDPYLHFFAYDPQYFLNELDLLGQAAAT
 QLARNISNSGAMQTYFAGEILPGDNLAYDADLSAWTEYIPYHFRPNYHDVGTCSMPKEMGSVVDNAARV
 YGVRGLRVIDGSIPPTQMSHVM TVFYAMALKISDAILEDYASMQ



Amplifikace genu GO2 z gDNA *Aspergillus niger*

Materiál:

Elizyme HS ROBUST MIX RED (Elisabeth Pharmacon)

Primery (viz. příložený Datasheet):

N-HIS,T7-GO2-HIS (T7-GO2-HIS_FW/GO2-HIS_REV)

N-HIS,T7-GO2 (T7-GO2-HIS_FW/GO2_REV)



N-HIS-GO2 (HIS-GO2_FW/GO2_REV)

Složení PCR reakční směsi (50 ul)

HS ROBUST MIX RED	25.0 ul
FW primer (10 uM)	2.0 ul
REV primer (10 uM)	2.0 ul
PCR voda	16.0 ul
gDNA	5.0 ul

Teplotní protokol:

95°C	2:30	1 x
95°C	10s	
60°C	20s	40 x
72°C	50s	
72°C	5 min	1 x



Přečištění PCR produktu

Materiál: Monarch® PCR & DNA Cleanup Kit (5 µg) (https://international.neb.com/-/media/nebus/files/protocols/t1030_quick_protocol_card_monarch_pcrdna_cleanup.pdf?rev=a562d8f1f21741b0ac6d59ea9020cef3&hash=E4845F80E8B2BE925A09DAD6562A41EC)

Restrikční analýza

Materiál:

Restriktázy

*Nde*I (<https://international.neb.com/products/r0111-ndei#Product%20Information>)

*Xho*I (<https://international.neb.com/products/r0146-xhoi#Product%20Information>)

*Bam*HI-HF (<https://international.neb.com/products/r3136-bamhi-hf#Product%20Information>)

Vektor pET28a

Složení směsi RA (50 ul)

CutSmart buffer (10x)	5.0 ul
PCR produkt/vektor	15 ul
PCR voda	28 ul
NcoI-HF/ BamHI-HF	1.0 ul (10U)
XhoI	1.0 ul (10U)

Teplotní protokol: Inkubace při 37°C po dobu 3 hodin

Izolace naštěpeného PCR produktu a plasmidu z agarózového gelu

Materiál: Monarch® DNA Gel Extraction Kit (https://international.neb.com/-/media/nebus/files/protocols/t1020_quick_protocol_card_monarch_dna_gel_extraction.pdf?rev=09308c01500f43c6a8589b01845765d9&hash=0FD4350A9DE2AEAC84913DB2A00140AF)

Ligace PCR produktu do vektoru pET28a

Materiál: Instant Sticky-end Ligase Master Mix

(<https://international.neb.com/protocols/2012/08/27/protocol-transfer-master-mix-to-ice-prior-to-reaction-set-up-mix-tube-by-finger-flicking-before-u>)

Protokol:

1. Transfer master mix to ice prior to reaction set up. Mix tube by finger flicking before use.
2. Combine 100 ng of vector pET28a with a 3-fold molar excess of insert and adjust volume to 5 µl with dH₂O.
3. Add 5 µl of Instant Sticky-end Ligase Master Mix, mix thoroughly by pipetting up and down 7-10 times, and place on ice. The sample is now ready to be used for transformation.

Transformace vektoru po ligaci do NEB 5-alpha E. coli

Materiál:

NEB 5-alpha Competent E. coli (High Efficiency) (<https://international.neb.com/products/c2987-neb-5-alpha-competent-e-coli-high-efficiency#Product%20Information>)

Protokol:

1. Thaw a tube of NEB 5-alpha Competent E. coli cells on ice for 10 minutes.
2. Add 2 µl of plasmid DNA to the cell mixture. Carefully flick the tube 4-5 times to mix cells and DNA. Do not vortex.
3. Place the mixture on ice for 30 minutes. Do not mix.
4. Heat shock at exactly 42°C for exactly 30 seconds. Do not mix.
5. Place on ice for 5 minutes. Do not mix.
6. Pipette 950 µl of room temperature SOC into the mixture.
7. Place at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.
8. Warm selection plates to 37°C.
9. Spread 50-100 µl of each dilution onto a selection plate and incubate overnight at 37°C.