

## 1.a. Potentiometric determination of dissociation constant of weak acid

---



The dissociation constant of weak  $HA$  is

$$K_A = \frac{a_{H^+} \cdot a_{A^-}}{a_{HA}} \quad (1.1.)$$

where  $a_{HA}$ ,  $a_{H^+}$  and  $a_{A^-}$  are activity of undissociated weak acid and activity of dissociation products respectively. By the mathematical rearrangement of the equation (1.1.), we obtain the Henderson-Hasselbalch equation:

$$-\log K_A = pK_A = pH + \log \frac{a_{HA}}{a_{A^-}} \quad (1.2.)$$

Activities can be replaced with molarities ( $a_{HA} = [HA]$ ,  $a_{A^-} = [A^-]$ ) in the case of dilute solutions below approx.  $10^{-2} M$  in this case.

The dissociation constant  $pK_a$  can be determined by titration of weak acid using a strong base or acid. Measurement of  $pH$  changes during titration allows evaluating the molarity of the dissociated and undissociated form of weak acid. Titration may be done in two ways, depending on the  $pK_a$  value. The way depends whether  $pK_a$  of weak acid is below or above value 7.

### TITRATION OF VERY WEAK ACID WITH STRONG HYDROXIDE.

Boric acid ( $HA$ ) represents a very weak acid  $pK_a > 7$ . The first degree dissociation constant  $pK_a$  can be determined from the titration results of boric acid with strong potassium hydroxide (optionally sodium hydroxide).

Electroneutrality and mass conservation law of weak acid are respected at any addition of strong  $KOH$  hydroxide to the weak  $HA$  solution:

$$[A^-] + [OH^-] = [K^+] + [H^+] \quad (1.3.)$$

$$c_{HA} = [HA] + [A^-] \quad (1.4.)$$

where  $c_{HA}$  means analytical concentration of  $HA$ . Potassium cation remains dissociated during titration and its concentration  $[K^+]$  is equal to the total analytical concentration of strong hydroxide  $c_{KOH}$ . Then the eqn (1.3.) gives:

$$[A^-] = c_{KOH} + [H^+] - [OH^-] \quad (1.5.)$$

We obtain for the undissociated acid concentration  $[HA]$  following if we insert the relation (1.5.) in the equation (1.4.) and then mathematically rearrange the result:

$$[HA] = c_{HA} - c_{KOH} - [H^+] + [OH^-] \quad (1.6.)$$

The obtained concentrations  $[A^-]$  and  $[HA]$  can be used in the Henderson-Hasselbalch equation (1.2.).

### TITRATION OF SALT OF WEAK ACID WITH STRONG ACID.

Acetic acid ( $HA$ ) represents a weak acid  $pK_a < 7$ . The dissociation constant  $pK_a$  can be determined from the titration results of salt  $KA$  (where  $K$  means sodium or potassium) of weak acid with strong acid  $HX$  (for example with hydrochloric acid  $HCl$ ).

The mixed solution of  $KA$  and  $HX$  keep electroneutrality and the mass conservation law of salt of weak acid:

$$[A^-] + [OH^-] + [X^-] = [K^+] + [H^+] \quad (1.7.)$$

$$c_{KA} = [A^-] + [HA] \quad (1.8.)$$

The salt  $KA$  is completely dissociated, so the concentration  $[K^+]$  during the titration is always equal to the total analytical salt concentration  $c_{KA}$ . For the same reason, the strong acid anion concentration  $[X^-]$  is equal to the total analytical concentration  $c_{HX}$  of strong acid.

If we replace the concentrations  $[K^+]$  and  $[X^-]$  in eqn (1.7.) by analytical concentrations, we get:

$$[A^-] = c_{KA} + [H^+] - [OH^-] - c_{HX} \quad (1.9.)$$

We can insert this relation in eqn (1.8.), mathematically rearrange and to obtain:

$$[HA] = c_{HCl} + [OH^-] - [H^+] \quad (1.10.)$$

The concentrations  $[A^-]$  and  $[HA]$  we use in the Henderson-Hasselbalch eqn (1.2.).



**TASK:** Determine  $pK_a$  of acetic acid and boric acid to the first stage. Compare the result with the values from literature.



**LABORATORY AIDS AND CHEMICALS:** automatic titrator equipped by  $pH$  detection, weighing boat, 2 volumetric flasks (50ml), 0,1M KOH or NaOH, 0,1M HCl, boric acid, sodium acetate, titration beaker (100ml).



**INSTRUCTIONS:** Get study how to operate the automatic titrator and how to measure the  $pH$ , switch titrator on and prepare the device for measurement. Check  $pH$  measurements with calibration solutions (use buffers).

- Weigh  $5 \cdot 10^{-4}$  moles of boric acid (or sodium acetate) with maximum accuracy. Transfer the weight to a 50 ml volumetric flask and dissolve in distilled water up to volume marking. Pour the solution into a dry titration beaker.
- Titrate with solution agent 0,1M KOH (or 0,1M HCl) by 0.5ml dose to total dose 5ml. Record the  $pH$  after stabilization.



**DATA ANALYSIS:** The ratio of the concentration of the dissociated and non-dissociated form  $[HA]/[A^-]$  of acid  $HA$  changes after each dose. The concentrations  $[H^+]$  and  $[OH^-]$  can be omitted in the data evaluation if they are significantly lower than other concentrations in eqns (1.5.), (1.6.), (1.9.) or (1.10.).



**REPORT:** For both very weak and weak acid: weight in mg, **TABULE 1:** for zero and each consequent titration dose: volume of added agent, total volume, analytical concentration of boric acid (or sodium acetate), experimental  $pH$ ,  $[H^+]$ ,  $[OH^-]$ ,  $[A^-]$ , and  $[HA]$  concentrations, expression  $\log[HA]/[A^-]$ ,  $pK_A$  calculated using eqn (1.2.). **NEXT:** mean value of  $pK_A$  and confidence interval using Student's t-distribution.