

Amplicon purification

Using AMPure beads to purify and size-select PCR amplicons

Prep:

- 01 Bring AMPure beads to **RT** before starting the purification (30 minutes).
- 02 **Prepare fresh 80% EtOH.** You will need approximately 0.5 mL per each sample.
- 03 **Vortex** the **AMPure** beads for at least **60 seconds**.
- 04 Combine the PCR reaction and the AMPure beads. Per 25 μL of PCR reaction you need to use **20 μL** of **AMPure** beads. **Mix properly** by pipetting or vortexing.
- 05 Incubate at RT for **5 minutes**. *Off the magnet.*
- 06 Put the samples on the **magnet**. Incubate for **2 minutes**.
- 07 Keep the samples on the magnet and remove the supernatant. Be careful not to damage the pellet.
- 08 Add **200 μL** of **80% EtOH**. Incubate for **30 seconds**. **Remove** the ethanol/supernatant.
- 09 **Repeat the previous step:** *Add 200 μL of 80% EtOH, wait 30 seconds, remove the EtOH.*
- 10 Close the tubes/put foil on the plate. Do quick spin to collect the remaining ethanol.
- 11 Put the samples back on the magnet. Wait a minute. Remove the residual ethanol with a 10 μL pipette. Be careful not to damage the pellet.
*The goal is to remove as much ethanol as possible – but **the pellet cannot dry up!***
- 12 Remove the samples from the magnet.
- 13 Add **32 μL** of **Elution** buffer.
- 14 Close the tubes/put foil on the plate. Mix the samples properly. Vortex. Quick spin.
- 15 Incubate at RT for **5 minutes**. *Off the magnet.*
- 16 Put the samples on the **magnet**. Incubate for **2 minutes**.
- 17 Transfer **30 μL** of the supernatant into a clean tube/plate. Be careful not to damage the pellet. *You can use the 10 μL three times to minimise the bead contamination.*