

# 1. Acid-base Equilibria

## 1.a. Photochemical determination of dissociation constant of acid-base indicator



The 3',3'',5',5''-Tetrabromo-m-cresolsulfonephthalein (bromocresol green) acid-base indicator behaves as a reversible system whose acidic form (yellow,  $HB^-$ ) changes into a basic form (blue,  $B^{2-}$ ) at pH range 3.8–5.4. The concentration of both forms of the indicator can be determined by the photometric method.

The univalent anion of the indicator dissociates according to the chemical equation:



*yellow solution*    *blue solution*

The thermodynamic equilibrium constant of the dissociation to the second degree is given by:

$$K_A = \frac{a_{H_3O^+} \cdot a_{B^{2-}}}{a_{HB^-}} \quad (1.2.)$$

where  $a_i$  ( $i = H_3O^+, B^{2-}, HB^-$ ) are the activities of the ions. The relationship between thermodynamic dissociation constant  $K_A$  and dissociation constant  $K'_A$  obtained from concentrations is:

$$K'_A = \frac{[H_3O^+][B^{2-}]}{[HB^-]} = K_A \frac{\gamma_{HB^-}}{\gamma_{H_3O^+} \cdot \gamma_{B^{2-}}} \quad (1.3.)$$

where  $\gamma_i$  are the activity coefficients of the ions. After mathematical rearrangement, we get:

$$pK'_A = pH - \log \frac{[B^{2-}]}{[HB^-]} \quad (1.4.)$$

The activity coefficients of ions can be obtained by use extended Debye-Hückel law (DHL). The activity  $\gamma_{B^{2-}}$  is given in aqueous solution at 25°C by expression:

$$\log \gamma_{B^{2-}} = -\frac{A \cdot (z_{B^{2-}})^2 \cdot \sqrt{I}}{1 + B \cdot r_{B^{2-}} \cdot \sqrt{I}} = -\frac{2,034\sqrt{I}}{1 + 2,30\sqrt{I}} \quad (1.5.)$$

where  $A = 0,5085$ ,  $B = 0,3281$ ,  $r_{B^{2-}} = 0,7 \text{ \AA}$  is effective diameter of the ion  $B^{2-}$  in Ångström. The ionic strength  $I$  is (at low concentrations) given by:

$$I = \frac{1}{2} \sum_{i=1}^k c_i z_i^2 \quad (1.6.)$$

where  $z_i$  are charge numbers of all ions  $i$  in the solution,  $c_i$  are their molarities.

The activity coefficients  $\gamma_{H_3O^+}$  and  $\gamma_{HB^-}$  are equal according DHL thus relationship between constants  $K_A$  and  $K'_A$  can be simplified to:

$$K_A = K'_A \gamma_{B^{2-}} \quad \text{ie: } pK^A = pK'^A - \log(\gamma_{B^{2-}}) \quad (1.7.)$$

and together with eqn (1.5.) it gives:

$$pK_A = pK'_A + \frac{2,04\sqrt{I}}{1 + 2,3\sqrt{I}} \quad (1.8.)$$

The thermodynamic equilibrium constant of the dissociation  $K_A$  can be calculated using eqn (1.8.) or it can be graphically evaluated more precisely from an experiment at different ionic strength.



**TASK:** Evaluate thermodynamic equilibrium constant of dissociation  $K_A$  of bromocresol green to the second degree at ionic strength  $0,1M$ .



**LABORATORY AIDS AND CHEMICALS:** UV/VIS spectrophotometer (minimum range  $350-720\text{ nm}$ ), 2 cuvettes, 2 volumetric flasks ( $50\text{ ml}$ ), 1 volumetric flask ( $250\text{ ml}$ ), 3 volumetric pipettes ( $1, 5, 25\text{ ml}$ ), 1 scale pipette ( $10\text{ ml}$ ),  $1,5 \cdot 10^{-4}M$  stock solution of bromocresol green (CAS No: 76-60-8),  $0,2M\text{ CH}_3\text{COONa}$ ,  $1M\text{ CH}_3\text{COOH}$ ,  $1M\text{ KCl}$ , and  $3M\text{ HCl}$ .




**INSTRUCTIONS:**

**Preparation of solutions I and II.** Prepare  $50\text{ ml}$  (use flask of same volume) of solution I of concentration  $1,5 \cdot 10^{-5} M$  bromocresol green (BG) inside  $0,01 M\text{ CH}_3\text{COONa}$  at ionic strength  $I=0,1M$  from stock solutions. Set the ionic strength to the desired value with a pre-calculated volume of  $1M\text{ KCl}$ . Prepare  $50\text{ ml}$  of solution II of concentration  $1,5 \cdot 10^{-5} M$  bromocresol green (BG) inside  $0,01 M\text{ CH}_3\text{COONa}$  at ionic strength  $I=0,1M$  from stock solutions.

**Measuring spectra of indicator at different pH.** Pour whole solution I into larger flask ( $250\text{ ml}$ ). Take a sample of solution I, place it in a quartz cuvette and measure the entire UV / VIS spectrum. Determine the wavelength at which the solution has a maximum absorbance  $A_2$  (see **FIGURE 1**). Return the content of the cuvette to the flask with the original solution I. Add  $1\text{ ml}$  of solution II to the flask and mix. The pH of solution is changed. Repeat sampling, spectrum measurement, sample return and addition of  $1\text{ ml}$  of solution II. Repeat this procedure 6 times. The last addition is done with  $1\text{ ml}$  of  $3M\text{ HCl}$ .

The solution containing the equimolar ratio  $\text{CH}_3\text{COONa}$  and  $\text{CH}_3\text{COOH}$  is green in color and has two maximas (see **FIGURE 1**).



**FIGURE 1:** Evaluation of spectra of acid-base indicator of bromocresol green obtained at different pH.



**DATA ANALYSIS:** The ratio of the concentrations of the basic and acid forms of the indicator is equal to the absorbance ratio at the adsorption maximum (compare **FIGURE 1**):

$$\frac{[B^{2-}]}{[HB^-]} = \frac{A_i - A_1}{A_2 - A_i} \quad (1.9.)$$

where  $A_2$  is absorbance of  $B^{2-}$  anion if the anion  $HB^-$  is not present (ie in a far basic environment).  $A_1$  is absorbance of  $HB^-$  anion if the anion  $B^{2-}$  is not present (ie in a far acidic environment),  $A_i$  is absorbance of  $B^{2-}$  anion at a general  $pH$  when both anions  $B^{2-}$ ,  $HB^-$  coexist in the solution.

The  $pH$  of the solutions to be monitored is determined by the concentration of the majority of the solution components, which are acetic acid and sodium acetate. They form conjugated acid-base buffer. The  $pH$  is given by Henderson-Hasselbalch eqn:

$$pH = pK^{HAc} + \log \frac{c^{NaAc}}{c^{HAc}} \quad (1.10.)$$

where  $pK^{HAc} = 4,76$  is the negative logarithm of the dissociation constant of acetic acid.  $c^{NaAc}$  and  $c^{HAc}$  are analytical concentrations of sodium acetate and acetic acid.



**REPORT: TABLE 1:** The volumes of the stock solutions used to prepare solutions I and II. The detailed calculation of the ionic strength. **Common graph 1:** UV/VIS spectra for all sample solutions. **Next:** wavelength of absorption maxima of  $B^{2-}$  and  $HB^-$ , value  $A_2$  and  $A_1$  (**FIGURE 1**). **Table 2:** for each sampling: addition of solution II, experimental absorbance  $A_i$ , calculated ratio  $(A_i - A_1)/(A_2 - A_1)$  (use eqn (1.9.)),  $\log[(A_i - A_1)/(A_2 - A_1)]$ ,  $c^{NaAc}$  and  $c^{HAc}$ ,  $pH$  value calculated using eqn (1.10.) and  $pK^{HAc}$  from literature.  $pK'_A$  (eqn (1.4.)),  $pK_A$  (eqn (1.8.)). **Next:** The mean value  $pK_A$  and its confidence interval according to the Student's t-distribution.