

**Not So Short Introduction To
TopSpin**

Spring 2021

Karel Kubíček

with the help of Pavel Veverka, Magdaléna Krejčíková,

and Radovan Fiala

Disclaimer:

This presentation is intended as a basic introduction to TopSpin 3.2 for those who are not familiar with it and the content wishes to cover initial steps to handle the software, sample and spectrometer properly, be able to acquire simple 1/2D NMR spectra and some hints for processing.

If imprecise or unclear formulations encountered, talk to the authors or to more experienced spectroscopist or facility staff.

Everybody is encouraged to study literature and discuss with more experienced colleagues to clarify all below written and become autonomous NMR operator.

Before the TopSpin runs OR some useful UNIX commands (in **blue**) and results (in **red**)

[karelk@nmrcf ~]\$ **whoami** - this command displays who's logged in

karelk

[karelk@nmrcf ~]\$ **hostname -a** - displays full name of the computer you're working on

nmrcf.ncbr.muni.cz

[karelk@nmrcf ~]\$ **quota** - displays space available for the current user

Disk quotas for user karelk (uid 16329):

| Filesystem | blocks | quota | limit | grace | files | quota | limit | grace |
|---------------------|--------|--------|--------|-------|-------|-------|-------|-------|
| /dev/nmrspect1/home | 25964 | 500000 | 600000 | | 1146 | 0 | 0 | |

[karelk@nmrcf ~]\$ **quota -g lrpi** - displays space available for the nmr data

Disk quotas for group lrpi (gid 2005):

| Filesystem | blocks | quota | limit | grace | files | quota | limit | grace |
|----------------------|----------|-----------|-----------|-------|-------|-------|-------|-------|
| /dev/nmrspect1/data1 | 53431604 | 100000000 | 110000000 | | 54667 | 0 | 0 | |

[karelk@nmrcf ~]\$ **pwd** - where are you working (print working directory)

/home/karelk

[karelk@nmrcf ~]\$ **cd /d1/data/karelk/nmr/** - change directory to /d1 ...

[karelk@nmrcf nmr]\$ **pwd**

/d1/data/karelk/nmr

This is one of the key directories for your NMR experiments, generally it is as follows:

/<disk>/data/<user>/nmr/<dataset> - the items in <> and italic (can) change depending on spectrometer and user. Datasets are of any name. Think of a good ordering system at the beginning of your NMR experience ...

/<disk>/data/<user>/nmr/<dataset>

```
[karelk@nmrcf nmr]$ cd 130118_700b_Dataset/ - change directory to a Dataset
```

```
[karelk@nmrcf 130118_700b_Dataset]$ ls - list directory contents
```

```
1 10005 2 3 4 5
```

```
[karelk@nmrcf 130118_700b_Dataset]$ grep PUL */acqu - print lines matching a pattern, e.g.
```

PUL – it is useful to see what kind of spectra were measured within the dataset, the output is however not sorted

```
10005/acqu:##$PULPROG= <c_con_iasq>
```

```
1/acqu:##$PULPROG= <zgpr>
```

```
2/acqu:##$PULPROG= <hsqcetf3gpsi2>
```

```
3/acqu:##$PULPROG= <c_caconcaco2_ct_nove>
```

```
4/acqu:##$PULPROG= <hccconhgp3d3>
```

```
5/acqu:##$PULPROG= <c_con_iasq>
```

```
[karelk@nmrcf 130118_700b_Dataset]$ grep PUL */acqu | sort -n - here it is sorted by  
number
```

```
1/acqu:##$PULPROG= <zgpr>
```

```
2/acqu:##$PULPROG= <hsqcetf3gpsi2>
```

```
3/acqu:##$PULPROG= <c_caconcaco2_ct_nove>
```

```
4/acqu:##$PULPROG= <hccconhgp3d3>
```

```
5/acqu:##$PULPROG= <c_con_iasq>
```

```
10005/acqu:##$PULPROG= <c_con_iasq>
```

/<disk>/data/<user>/nmr/<dataset>

```
[karelk@nmrcf 130118_700b_Dataset]$ cd 1 - change directory to Experiment number 1 (ExpNo 1)
```

```
[karelk@nmrcf 1]$ ls
```

```
acqu  acqu  audita.txt  fid  format.temp  pdata  pulseprogram  scon2  
shimvalues  specpar  uxnmr.info  uxnmr.par - Note that there are several files and directory  
pdata.
```

acqu – contains acquisition parameters prior measurement

acqu – contains info about the parameters as they were at the end of measurement

audita.txt – info about who/where/when ran the measurement

fid – raw one dimensional data

pulseprogram – compiled pulseprogram

pdata – directory containing processed data in different processing number directories (ProcNo, mostly 1)

```
[karelk@nmrcf 1]$ cd pdata
```

```
[karelk@nmrcf pdata]$ ls
```

```
1
```

```
[karelk@nmrcf pdata]$ cd 1/
```

```
[karelk@nmrcf 1]$ ls - 1r and 1i are the real and imaginary part of your processed 1D spectrum
```

```
1i  1r  assocs  auditp.txt  outd  proc  procs  thumb.png  title
```

```
[karelk@nmrcf 1]$ pwd
```

```
/d1/data/karelk/nmr/130118_700b_Dataset/1/pdata/1
```

/<disk>/data/<user>/nmr/<dataset>/<ExpNo>/pdata/<ProcNo>/ - where
ExpNo and ProcNo may be integer 1..9999 (may be more but never tested)

/<disk>/data/<user>/nmr/<dataset>

Let's move back to the Dataset directory and change directory to nD experiment, here ExpNo 2

```
[karelk@nmrcf 130118_700b_Dataset]$ pwd
```

```
/d1/data/karelk/nmr/130118_700b_Dataset
```

```
[karelk@nmrcf 130118_700b_Dataset]$ cd 2
```

```
[karelk@nmrcf 2]$ ls
```

```
acqu  acqu2  acqu2s  acqus  audita.txt  cpdprg3  format.temp  pdata  
pulseprogram  scon2  ser  shimvalues  specpar  uxnmr.info  uxnmr.par
```

Note that there is similar content as of 1D experiment with two important differences – **fid** is replaced by **ser** and as there are two dimensions, there are **acqu2** and **acqu2s** files.

```
[karelk@nmrcf 2]$ cd pdata/1/
```

```
[karelk@nmrcf 1]$ ls
```

– in case of a 2D spectra, the processed data are in 2rr. Practically speaking one can ignore the imaginary files once the spectrum is properly processed and phased. Keeping just the 2rr file saves disk space.

```
2ii  2ir  2ri  2rr  assoc  auditp.txt  clevels  outd  proc  proc2  proc2s  
procs  thumb.png  title
```

```
[karelk@nmrcf 1]$ pwd
```

```
/d1/data/karelk/nmr/130118_700b_Dataset/2/pdata/1
```

/<disk>/data/<user>/nmr/<dataset>/<ExpNo>/pdata/<ProcNo>/

/<disk>/data/<user>/nmr/<dataset>

Let's move back to the Dataset directory and change directory to another nD experiment, this time ExpNo 4 and the pdata to see ProcNos different from 1

```
[karelk@nmrcf 130118_700b_Dataset]$ pwd
/d1/data/karelk/nmr/130118_700b_Dataset
[karelk@nmrcf 130118_700b_Dataset]$ cd 4/pdata
[karelk@nmrcf pdata]$ ls
1 13 23 998 999
```

This is typical content for a 3D dataset where in ExpNo 1 there is the processed 3D, as one has dimensions 1-2-3, ProcNo 13/23 are first planes of the spectral cube in 13/23 direction. 998 and 999 represent that any ProcNo is allowed.


/<disk>/data/<user>/nmr/<dataset>/<ExpNo>/pdata/<ProcNo>/

NMR Acquisition


TopSpin

Starting TopSpin – A) Application menu or B) command line


A) Application Menu



B) Terminal / command line



B) Terminal / command line



```
karelk@nmr600:~  
File Edit View Search Terminal Help  
[karelk@nmr600 ~]$ topspin  
TopSpin 4.0.9 - Copyright (C) 2020 Bruker BioSpin GmbH  
Installation directory: "/opt/topspin4.0"  
Waiting for CodeMeter license  
█
```



```
karelk@nmr600:~  
File Edit View Search Terminal Help  
[karelk@nmr600 ~]$ topspin  
TopSpin 4.0.9 - Copyright (C) 2020 Bruker BioSpin GmbH  
Installation directory: "/opt/topspin4.0"  
Waiting for CodeMeter license  
valid until 2035-11-14 14:00:08 +0100  
█
```

topspin 

NMR With Ease


Connecting to command service...

The Next Generation in NMR Software

© 2020 Bruker BioSpin

Topspin 4.X – layout description


- 1 – Menu bar
- 2 – Spectra navigator
- 3 – Toolbar
- 4 – Browser / history
- 5 – Data/window area
- 6 – Command line
- 7– Status panel



Right mouse click in the status bar area toggles it on/off


Topspin 4.X – layout description

- 1 – Menu bar
- 2 – Spectra navigator
- 3 – Toolbar
- 4 – Browser / history
- 5 – Data/window area
- 6 – Command line
- 7 – Status panel



Topspin 4.X – layout description


Browser / history



Search: Find

• /d1

No structure available.




Search: Find

- karelk
 - 1D1H
 - 1D2H
 - 030409_CID
 - 081127_SS_Dach_288K
 - 090113_SS_Dach_278K
 - 090126_SS_Dach_278K
 - 090402_CID_titr
 - 090421_SS_Dach_Imino278K
 - 090609_SS_Dach_Imino288K
 - 090703_CID_Noesy13C
 - 091203_SS_Dach_Imino278K
 - 100304_CID_Titr
 - 100309_CID-CTD_HCCH-TOCSY
 - 100312_CID-CTD_HNHa_CbHarom
 - 100331_PAZ
 - 100413_AgoPAZ

Status Line History

| Time | Message |
|----------|--|
| 17:41:49 | Connection: searching TOPSPIN comm... |
| 17:41:49 | Connection: TOPSPIN system found... |
| 17:41:49 | Login to CPR in progress ... |
| 17:41:49 | Logon accepted. |
| 17:41:50 | Login to CPR succeeded. |
| 17:41:50 | dataserver: started |
| 17:41:50 | dataserver: in progress |
| 17:41:50 | dataserver: Server 'corbaloc:iiop:147.2... |
| 17:41:50 | Connecting to restartserver |
| 17:41:51 | Connecting to restartserver |
| 17:43:25 | 2021-03-08 11:29:46.779: Spectromet... |
| 17:43:25 | 2021-03-08 11:29:46.780: Device Blad... |
| 17:43:25 | 2021-03-08 11:29:46.793: Power outl... |
| 17:43:25 | 2021-03-08 11:31:12.856: Spectromet... |
| 17:44:00 | Starting BSMS display... |





Search: Find

- karelk
 - 1D1H
 - 1D2H
 - 030409_CID
 - 081127_SS_Dach_288K
 - 090113_SS_Dach_278K
 - 090126_SS_Dach_278K
 - 090402_CID_titr
 - 090421_SS_Dach_Imino278K
 - 090609_SS_Dach_Imino288K
 - 090703_CID_Noesy13C
 - 091203_SS_Dach_Imino278K
 - 100304_CID_Titr
 - 100309_CID-CTD_HCCH-TOCSY
 - 100312_CID-CTD_HNHa_CbHarom
 - 100331_PAZ
 - 100413_AgoPAZ

Command Line History

| Time | Message |
|----------|----------|
| 17:44:00 | bsmsdisp |






Spectrum display properties.

Toggle the visibility of integrals, peaks, and other spectra components.

Note: Any of these components can only be displayed when available, e.g. when peaks were picked before, when integrals were calculated, when a structure mol file is available, etc.



Notes to TopSpin

When you open Topspin for the first time under your own account, you need to set up a few things

- Set the path to data
- Create directory for your data
- Create a dataset in your data directory
- Open and customize the status panel
- Disable zg safety
- Other customizations available – colors, # of command lines, lock grid lines etc.

The preferences are stored in your home directory, e.g.

`/home/fiala`

In the directory `.topspin1`

May contain some junk, the essential file is




`/home/fiala/.topspin1/prop/globals.prop`

If you want to use someone else's setup, copy this directory/file from his/her home directory to your home directory before starting TopSpin



- Prior any manipulation with the spectrometer or TopSpin, check the status panel first – namely, what is the acquisition status, temperature and what is the lock signal status.

| | | | | | | | | | | | |
|---------------------|-------------------|-------------------------|-----------|------|--------|-----------------------|--------|-------------------------------|-----------------------|------------------------------------|------|
| Spectrometer Status | Amplifier Control | Acquisition information | Fid Flash | Lock | Sample | Shim Coil Temperature | POWCHK | Probe Temperature | Temperature Channel 1 | Spooler | B |
| On ☑ | | no acquisition running | | | | 301 K | ✓ | 298.0 K On ☑ Reg. State: ☑ | 298.0 K | queued: 0 delayed: 0 cron: 0 | Autc |

- In case there is a sample of your colleague still in the magnet, take it out from the magnet (*vide infra*) and take best care about the sample.





- To eject the sample from the magnet – assure yourself there is no acquisition running, the lock is off (no straight line is sweeping the lock window but dispersive sinusoidal curves can be seen in the lower area of the lock window – see the figure on the right) and there is no mechanical obstacle that could prevent sample ejection. Then either type **ej** on the cl or use the BSMS window to manipulate with the lift.



Ejecting the sample

- In the Main panel of the BSMS window (Bruker Smart Matching /Shimming system), click Sample - Lift for turning on the airflow (button turns green)
- Sample mail will be set in action and brings the sample down
- Spectrometers that have no sample mail will announce the sample ejection by increased airflow and sample „dancing“ on the top of the magnet. In this case, after removing the sample from the magnet, do not forget to switch off the airflow by clicking the Lift button once again (turns grey)



Temperature control

- Once there is no sample in the magnet, set the temperature you want for your sample/measurement – this is done by opening the Temperature control window. This can be done by *i)* double click the Sample Temperature window in the status panel, *ii)* type **edte** in the command line (cl). This will open new window called **T**. Set the temperature and in the monitoring tab check first four check-boxes and control the progress of temperature. As soon as the temperature is stabilized, you may insert your sample into the magnet.

Temperature Monitoring Record Correction Self tune Configuration Log Help

On Off VTU State: On

| Channel | Regulation State | Stability | Current Temperature | Target Temperature | Heater Power |
|--|---|--|------------------------------------|------------------------------------|-----------------------------------|
| 1 5 mm CPQCI 1H-31P/13C/15N/D Z-... | <input checked="" type="checkbox"/> Steady | <input checked="" type="checkbox"/> Stability Lost | 293.2 K | 293.2 K (233.0 K...353.0 K) | 11.9 % (max. 35.0 % of 42.9 W) |
| | | | | <input type="button" value="Set"/> | |
| | State | Gas Flow | Target Gas Flow | Standby Gas Flow | |
| Probe Gas | <input checked="" type="checkbox"/> Steady | 400 lph | 400 lph | 200 lph | |
| | | | <input type="button" value="Set"/> | <input type="button" value="Set"/> | |
| Accessory Channel | State | Current Power | Target Power | | |
| 1 (Chiller) BCU | <input checked="" type="checkbox"/> Connected | Medium | Medium | | |
| | | | <input type="button" value="Set"/> | | |

Stability Lost 293.2 K (233.0 K...353.0 K)

Gas Flow Target Gas Flow Standby Gas Flow

400 lph 400 lph 200 lph

Set target temperature

Please enter the new probe target temperature.

Target temperature [K]:

Temperature Control Suite

Temperature **Monitoring** Record Correction Self tune Configuration Log Help

Configure



Temperature display properties: Show overview Auto scale

Power display properties: Show overview Auto scale

Gas flow display properties: Show overview Auto scale

General display properties: Update interval [s]: 1

| Channel | Current Temperature | Target Temperature | Current Power |
|--|---------------------------------------|---------------------------------------|---------------------------------------|
| 1 5 mm CPTCI 1H/19F-13C/15N/D Z-... | <input checked="" type="checkbox"/> - | <input checked="" type="checkbox"/> X | <input checked="" type="checkbox"/> - |




Sample/tube handling

- Inserting the tube into the spinner
 - Do not push the tube straight down the spinner, the tube may break down
 - Screw the tube slowly into the spinner and at the same time gently push the spinner upward the tube

Adjusting/centering the tube in the spinner

- Put the tube with the spinner into a depth gauge and screw the tube until the length of the sample is symmetrical around the middle line crossing the coil-representing boxes




Notes to temperature control

- Temperature is maintained by flowing cooled/heated gas around the sample.
- Sample temperature may differ from the set temperature (sensor not in the sample, calibration of the sensor itself) – calibration needed
- Calibrations based on the change of chemical shift with temperature
- Methanol (neat, 4%, 0.2%), 25°C and down
- Ethylenglycol, 30°C and up
- NMR thermometer – calibration substance added to buffer, ^2H spectra measured during the experiment run, special locking parameters must include the calibration data

Inserting the sample

- After the desired temperature is reached, in the Main panel of the BSMS window (Bruker Smart Magnet System), click (skip this in case of sample mail) Sample - Lift for turning on the airflow (button turns green)
- Place the sample either into the sample mail or on the top of the magnet bore when „maximum“ airflow is reached
- Turn off the airflow by clicking the Lift button again (goes gray again)
- For the magnets equipped with sample mail, place the spinner with the tube into the sample rack and close it
- Wait until the sample is positioned in the probe



Locking the magnetic field

- Locking on the reference nuclei (deuterium)
- Once the sample is in the magnet (check Sample status icon on status bar, status LED of the sample mail etc.), usually lock signal appears as red/green signal in dispersion mode. In case you do not see that line – it may indicate that previous user was using different deuterated solvent or there is no ^2H in your sample. Undocumented **lopo** command reads in the lock parameters without engaging the lock.
- Type **lock** on the cl

Note: you may not be able to lock if the field homogeneity is poor. It is a good idea to read in suitable shims before you try to lock (command **rsh**). If you do not have suitable shims, look for a recent sucrose (for water samples) or lineshape (organic solvents) files.

Locking the magnetic field

- Table with solvents will pop up
- Select the right solvent from the table and click the OK-button or double click the solvent of choice

The screenshot displays the Bruker software interface with the 'Solvents table' dialog box open. The dialog box contains a table with the following data:

| Solvent | Description |
|--------------------|--|
| Acetic | acetic acid-d4 |
| Acetone | acetone-d6 |
| C6D6 | benzene-d6 |
| CD2Cl2 | dichloromethane-d2 |
| CD3CN | acetonitrile-d3 |
| CD3CN_SPE | LC-SPE Solvent (Acetonitrile) |
| CD3OD_SPE | LC-SPE Solvent (Methanol-d4) |
| CDCl3 | chloroform-d |
| CH3CN+D2O | HPLC Solvent (Acetonitril/D2O) |
| CH3OH+D2O | HPLC Solvent (Methanol/D2O) |
| D2O | deuteriumoxide |
| D2O_salt | deuteriumoxide with salt |
| Dioxane | dioxane-d8 |
| DMF | N,N-dimethylformamide-d7 |
| DMSO | dimethylsulfoxide-d6 |
| EtOD | ethanol-d6 |
| H2O+D2O | 90%H2O and 10%D2O |
| H2O+D2O_salt | 90%H2O and 10%D2O with salt |
| HDMSO | 90%DMSO and 10%DMSO-d6 |
| Juice | fruit juice |
| MeOD | methanol-d4 |
| Plasma | blood plasma |
| Pyr | pyridine-d5 |
| T_H2O+D2O+Me4NCI | (CD3)4NCI in 90%H2O and 10%D2O, for NMR thermometer |
| T_H2O+D2O+NaAc | sodium acetate in 90%H2O and 10%D2O, for NMR thermometer |
| T_H2O+D2O+Pivalate | pivalate-d9 in 90% H2O and 10% D2O, for NMR thermometer |
| T_MeOD | methanol-d4, for NMR thermometer |
| TFE | trifluoroethanol-d3 |
| THF | tetrahydrofuran-d8 |
| Tol | toluene-d8 |
| Urine | urine |

The background interface shows the 'Acquire' and 'Process' tabs, a search bar, and a status bar at the bottom. The status bar includes the following information:

- Spectrometer Status: On
- Amplifier Control: no acquisition running
- Acquisition information: 171012_600_C9320B 2 1 /d1/data/karelk/nmr
- Fid Flash: Lock
- Sample: Shim coil temperature: 301 K
- POWCHK: On
- Probe Temperature: 298.0 K
- Temperature Channel 1: 298.0 K
- Spooler: queued: 0, delayed: 0, cron: 0

1D *2 3D /2

- Search: Find
- 151210_600_000001_p11_short
 - 151217_600_TOAC2xpY1CTD_check
 - 151217_600_clamp7_N98C_S102C_C181
 - 151222_600_Nterm_SH2
 - 151222_600_titr_TOAC_tSH2_sumo_out
 - 151227_600_titr_TOAC_tSH2_sumo_out
 - 160406_600_Spt6_pY_SingleAAtitration
 - 160414_600_aRRM_Nab3
 - 160415_600_C7995
 - 160419_600_aRRM_Nab3_titr
 - 161109_600_C9230F_Ubiq
 - 171012_600_C7995_Practical1
 - 171012_600_C9320B
 - 1 - zgpr - zgpr
 - 2 - **hsqcetf3gpsi - HSQCETF3GPSI**
 - C7995
 - CID_900
 - CID_950S

Main Lock/Level Shim Autoshim Service Log Help

| Lock | Phase | Power | Gain | Shim |
|------|-------|-------|------|------|
|------|-------|-------|------|------|

| | | | | |
|--------|-------|-------|------|--|
| LOCK | | | | |
| On-Off | Phase | Power | Gain | |

| | | | | |
|--------|------|---------|------|---------------|
| SAMPLE | | | | |
| LIFT | SPIN | Measure | Rate | Always Locked |

| | | | | |
|----------|--------------------------------|----------------|----------------|----------------|
| SHIM | | | | |
| Spin: | Z | Z ² | Z ³ | Z ⁴ |
| NonSpin: | X | XZ | | |
| | Y | YZ | | |
| | XY | | | |
| | X ² -Y ² | | | |

STD BY

Absolute Difference

Previous Actual Step

Reset

Stepsize

STD BY

Shut down BSMS

Shut down BSMS Shut down BSMS before switching BSMS off

Config

External

Sample: down missing up Shim coil temperature 301 K

| Time | Status Message |
|----------|---|
| 17:54:30 | xf2: in progress |
| 17:54:30 | xf2: finished |
| 17:54:31 | topshim: 1 set shims |
| 17:54:31 | topshim: 2 target 1, fit 2 |
| 17:54:32 | topshim: 3 acquire data |
| 17:54:33 | stored NS 4 into /root/topspin-Blade... |
| 17:54:33 | topshim: 3 2D Fourier transform |
| 17:54:33 | xf2: in progress |
| 17:54:33 | xf2: finished |
| 17:54:33 | topshim: 3 process field map |
| 17:54:33 | topshim: 1 sample size = 2.21 cm, p... |
| 17:54:33 | topshim: 3 stop acquisition |
| 17:54:34 | stored NS 4 into /root/topspin-Blade... |
| 17:54:34 | acquisition finished |
| 17:54:34 | zg finished |
| 17:54:34 | topshim: 2 finish |
| 17:54:33 | topshim: 3 lock status |
| 17:54:4 | topshim: completed successfully |

Lock Display


| | | | | | | | | | |
|---------------------|-------------------|-------------------------|-----------|------|--------|-----------------------|--------|-------------------|-------------|
| Spectrometer Status | Amplifier Control | Acquisition information | Fid Flash | Lock | Sample | Shim Coil Temperature | POWCHK | Probe Temperature | Temperature |
| On | | no acquisition running | | | | 301 K | ✓ | 298.0 K | 298.0 K |

Homogeneity of the field

Adjusting homogeneity of the magnetic field – as the more homogeneous magnetic field results in a narrower lock signal which results in a higher d.c. voltage, one aims for an optimum lock signal by adjusting various shim currents

Shimming options

- Manual (in the BSMS display via the Shim tab)
- Automatic
 - Gradient shimming – Topshim (additional parameters, 3D only works with water samples)
 - Simplex (slow, works based on lock level, shims to optimize must be specified in a file)
 - Hardware autoshim (Autoshim tab in BSMS), activate before starting long experiments .
- As there may be needed to adjust the lock phase, lock gain and other parameters, type **loopadj** on the cl (optimizes lock phase, lock gain, loop phase, ...) which will take care about all of it. In case there is sufficient lock signal, menu with three lines will appear to confirm everything went smoothly and was set up. Should there be low lock signal, error message will appear.




Selecting the nuclei

- Prior starting any set up, pulse checks or measurements, make sure the channel routing is correct and all the channels that are going to be measured are active.
- To do so, type **edasp** on the cl
- set the desired nuclei from the pop-up menus: typically channel 1 – ^1H , channel 2 - ^{13}C , channel 3 – ^{15}N
- Click on Default to set proper routing

for the channels and Save and Close

the settings.



Tunning and matching the probe

- For obtaining optimal signal-to-noise ratio, one needs to tune the probe with the sample inserted. It is done by adjusting two mutually interactive capacitors. One tunes the circuit to the desired resonance frequency (**tuning**) and the other matches the impedance (**matching**). Tuning/Matching can be achieved either manually directly on the probehead that is NOT equipped with ATM unit or automatically from a PC in case the probehead IS equipped with the ATM unit. Automatic Tuning/Matching Automatically (atma)
- Start with **atma** (AutoTuneMatchAuto) on the cl – typically, this will tune the probehead automatically without any human intervention
- Once atma is done, it is recommended to check the tuning manually
- Type **atmm** (Automatic Tuning/Matching Manual) on the cl – window with a black wobbling curve that needs to be ftuned to the minimum indicated by vertical red line is opened
- The precision of the displayed curve is driven by wobble sweep width. Set it prior wobbling by command **wbsw** on the cl to 2 MHz or click wbsw icon and set it there.
- N.B. Sometimes the atmm gets stuck when trying to change the sweep width. Using of following command and its option **atmm manwbsw** has proven to improve the wobble behavior. Again, create own macro or use **wkk**.


Tuning and matching the probe

- While tuning the probehead „manually “ with atmm, keep an eye on the lock-level signal. It may drop significantly during the tuning so once finished, run topshim once again. Typically one is tuning the channels starting at lowest frequency (typically, ^{15}N followed by ^{13}C and then ^1H) but in case of (cryo)-probeheads with a new design (lock signal is dropping while tuning/matching) are supposed to be tuned from highest frequencies down (i.e. ^1H -> ^{13}C -> ^{15}N). In any case, check the tuning curve of a given nuclei iteratively.
- Tune/match each channel starting from the nuclei with the lowest frequency (on the 850/950MHz spectrometers in reversed order)
- Adjust the displayed sweep width (parameter WBSW (wobble sweep width) in the upper right corner, *vide infra*)
- Get the minimum of the curve to required position by clicking the arrows that represent expected step of the tune/match
- Start with matching (get the minimum to the bottom), then tune (center the minimum on the reference line). As tune/match are interconnected, the procedure becomes iterative but optimum has to be always reached.
- For samples with high salt [$c(\text{NaCl}) \geq 250 \sim 300 \text{mM}$], use either shaped NMR-tube or 3mm tube. Both tubes require special rotors and extremely gentle handling!

SPECTRUM PROCPARS ACQPARS TITLE


171 change wobble sweep width [wbsw] a/karelk/nmr

zgpr



Creating a dataset

- Steps involving edasp and atma/atmm require an existing dataset with parameters. Users typically keep in their data-tree directorie Setup/Calibration or something similar. In those dataset one usually doesn't change more than power-levels, (de)activates nuclei but not much more ...
- If one wants to create a dataset, there are several ways to do it
 - Copy an existing dataset you or someone else used previously
 - Select a dataset that worked previously, ideally by clicking right button in the browser and open in a new window. Often, different set is selected in browser and different in the window. Always check which spectrum you are working with!!!
 - Read in a standard parameter set (rpar)



Creating a dataset

- Once being in the right dataset, type **edc** on cl and set the new destination of the experiment, experimental number, title, user ...
- In case dataset should be copied to the next ExpNo, use command **ixpno** which will copy the dataset to next ExpNo

Create New Dataset - new

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the Advanced field.

Dataset

NAME: 171
EXPNO: 3
Directory: /d1/data/karelk/nmr
 Open in new window

Parameters

Use current parameters
 Read parameterset [] Select
 Set solvent: H2O+D2O
Additional action
 No additional action
 Execute getprosol
 Keep parameters: P 1, O1, PLW 1 Change

Advanced

Number of datasets (receivers): 1
Title: zgpr

OK Cancel More Info... Help

Create New Dataset - new

Prepare for a new experiment by creating a new data set and initializing its NMR parameters according to the selected experiment type. For multi-receiver experiments several datasets are created. Please define the number of receivers in the Advanced field.

Dataset

NAME: 171012_600_C9320B
EXPNO: 3
Directory: /d1/data/karelk/nmr
 Open in new window

Parameters

Use current parameters
 Read parameterset: ZGPR Select
 Set solvent: H2O+D2O
Additional action
 No additional action
 Execute getprosol
 Keep parameters: P 1, O1, PLW 1 Change

Advanced

Number of datasets (receivers): 1
Title: zgpr

OK Cancel More Info... Help


- If creating a new dataset, type **rpar** on cl to read standard dataset parameters of a required pulse sequence
- Type **getprosol** on cl to set actual pulse calibration. Should calibration of any nuclei vary from the standard ones (typically ^1H), type **getprosol 1H 13.55 13W** which means that 90° ^1H pulse is $13.55\mu\text{s}@13\text{W}$ taken as reference.
- Modify the rest of the parameters in AcqPars in Data area (TD – number of points collected during FID, DS – number of dummy scans, NS – number of scans, SW – spectral width, O1P – carrier frequency, ...)

Calibration of 90° ¹H pulse

- Automatic calibration
 - Type **pulsecal** on cl – info-panel will pop-up with the length and power of 90° pulse. The parameters will be also automatically updated within the dataset
 - Pulsecal gives results reasonable enough to use them for nD experiments or to use them as a starting point for manual calibration
- Manual calibration
 - Type **p1** on cl, set the pulse length 1 μs and type **zqgfp**. This will run a macro that is composed of **zg** (zero go - start an experiment), **qsine** (multiply resulting FID with QSINE window function) and **fp** (perform Fourier transformation with phase correction). Alternatively **zgefz** does the same but instead of qsine exponential window function is used for FID apodization.
 - Water-line signal will appear at about ~4.7 ppm. Phase the line to pure absorptive line either manually through Process menu and manual phase adjusting. Phasing window can also invoked by typing **.ph**. If there are few signals in the spectrum automatic phase correction – **apk** command will work well. Check also symmetry of the signal and/or any abnormalities of the line-shape. Also measure the LWHH to get an estimate about the quality of the shimming.
 - Proton hard pulse (90° or $\pi/2$) is strongly dependent on the tube diameter, salt concentration, temperature and solvent used. The pulse-length can be as short as 8 μs but can reach ~17 μs for high salt samples. As the 90° pulse is usually determined by measuring the 360° since the pulse gives a minimum signal, type **p1** on cl and set it to four times the expected length of 90° ¹H hard pulse (i.e. 32 to 68 μs). From the knowledge of sine-function one can easily guess whether the signal is longer 360° or shorter.
 - Once minimum in 360° is achieved, take one fourth the length of the pulse that gave a zero signal to obtain length for the 90° pulse.
 - Samples in organic solvents or D₂O can be calibrated on 180° pulse (first pass through zero).

Calibration of 90° ¹H pulse

Another option how to calibrate pulse in to use command **paropt**, which will pop-up a window where one sets up initial value [of the pulse], increment, and number of increments. Result of such paropt will provide similar output (initial value: 1μs, increment:1μs, number of increments: 50).




Optimizing acquisition parameters

- Before starting any acquisition, it is wise to check whether all parameters are at their optimum or can be further improved. For the sake of optimization there is command **gs** (go scan) that will run an “infinite” loop enabling real-time manipulation with pulse-lengths, delays, receiver gain etc. Immediate impact on the FID is observed.
- Start with optimizing receiver gain – either by typing **rga** prior running **gs** or manually till ADC overflow disappears*
- The parameters can be adjusted by:
 - dragging the slider (not recommended)
 - directly typing desired value
 - clicking on the slider-scale (clicking will be affected by the sensitivity)


ADC range lines

*note: cryoprobes tend to give false ADC overflow warnings




Optimizing acquisition parameters

- Other parameters that could be set are the exact position of receiver (Offset-tab), phase correction (particularly in case of WATERGATE water suppression; CorPhase) and duration of calculated shaped pulses (Pulse-tab).
- Always check the Warning message and FIDAREA which one wants to reach as low as possible as the signal in biomolecular samples is dominated by water which is supposed to be efficiently suppressed.



Start of the experiment

- Duration of the experiment
 - Type **expt** (experimental time) – gives the experimental time, required space on disc for the current dataset
 - Type **multiexpt** – allows to estimate duration of multiple experiments that are in row (e.g. ExpNos 1, 2, 3, 4 ...). Returns the time and date, when the experiments are finished.
- Running the experiment
 - Type **zgqfp** – starts and processes 1D experiment (composite macro of zg+qsin+fp) – in case of low concentrated samples, one set number of scans as high as 128 or even 1024. As such experiment takes already a significant amount of time, it is wise to check whether something is being acquired. Either check the on-line Fourier transform in the acquisition window or type **tr** to transfer the so far acquired data from spectrometer to computer and process them with **efp** or **qsin+fp**.




- Type **zg** – starts the experiment without processing it
- Type **stop** – immediately stops the experiment
- Type **halt** – stops the experiment after the current scan

Start of the experiment

- In case of an nD experiment follow the same steps as in case of 1D, i.e. first set all parameters properly and check with gs-command that everything is correct and is not hurting sample/probe. Before typing zg, check the duration with expt or multiexpt to see when the experiment is going to finish.
- When the experiment is started with the zg-command, wait till the first FID is acquired. In the mean-time check that the temperature and lock are stable and not affected by the running experiment.
- Once first FID measured and stored, type **rser 1**, which will transfer the first FID to ~TEMP directory ExpNo 1, ProcNo 1. Process the FID with qsin+fp and phase it. This will give you a rough information about signal/noise ratio of your experiment – if there is no signal in the first row, there will be not much signal at the end of the experiment (exceptions are, e.g. HNCACB or DQF COSY, but these are not experiments for beginners:-).

Start of the experiment

- Should it be possible to compare the running experiment with previous dataset, store the current ProcNo to ProcNo 2 by typing **wrp 2 y** – write ProcNo to 2 and if there is already ProcNo 2, **y**, overwrite it. Then go to the previous/reference experiment, type again rser 1; qsin; fp, phase the spectrum. Type **.md** or click the multiple spectra icon and then type **rep 2** – read ProcNo 2. Of course one can do all the overlay with mouse only but the clicking may be more time consuming as the directory gets filled with experiments and ~TEMP remains at the top (scrolling up, double-clicks ..., keyboard is keyboard:-).



Quick reference:

- 1) bsmsdisp – open control panel
- 2) ej – eject sample
- 3) edte – open temperature control panel
- 4) ij – inject sample
- 5) lockdisp – display lock window
- 6) lock – select solvent for locking the magnet
- 7) topshim – automatic shimming
- 8) loopadj – adjust lock parameters
- 9) edasp – set up spectrometer routing
- 10) wbsw – set up wobble sweep width
- 11) atma – automatic tuning/matching
- 12) atmm – manual -----”-----
- 13) edc – copy dataset
- 14) iexpno – copy dataset to next ExpNo
- 15) rpar – read parameter set
- 16) getprosol – set up pulses according to prosol
- 17) pulsecal - calibrate ^1H 90° pulse
- 18) paropt – optimize parameter (e.g., p1 length)
- 19) rga – automatically set up receiver gain
- 20) gs – go scan – optimize acquisition parameters
- 21) expt – estimate duration of the experiment
- 22) multiexpt - ----”---- for multiple expts in row
- 23) zg – start acquisition
- 24) zgqfp – acquire 1D and apply qsin apodization
- 25) qsin – multiply FID with qsin window function
- 26) xfb – process 2D
- 27) xfb n – process 2D and remove 2ri, 2ir, 2ii
- 28) tr – transfer data 1D from spectrometer to PC
- 29) edmac *name* – create or edit macro *name*
- 30) wrp 2 y – write processed data to ProcNo 2 and if exists, overwrite
- 31) rser 1 – extract first FID of an nD experiment
- 32) .ph – start phase menu
- 33) .md – spectra overlay window
- 34) apk – automatic phase correction
- 35) show – show active processes
- 36) curplot – set up printer
- 37) print – print spectrum
- 38) acqu – switch to acquisition window
- 39) ii – if spectrometer doesn’t communicate
- 40) ii restart – if ii doesn’t help
- 41) stop – stop immediately acquisition, rough
- 42) halt – stop it smoothly, recommended
- 43) pulse – calculate pulse length based on 90° hard pulse parameters
- 44) calcpowlev – similar to pulse