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

Application/performance and technical specifications	Biacore T100 – Unmatched performance	Biacore A100 – Unmatched productivity	Biacore Flexchip – Array-based parallel kinetic profiling	Biacore X100 – Ready-to-run research system	Biacore 3000 – Interaction analysis with SPR-MALDI interface	Biacore C – Rapid and reliable protein quantification
						
Application						
Kinetic/affinity characterization	● ● ●	● ● ●	●	● ● ●	● ●	-
Kinetic/affinity screening and profiling	● ●	● ● ●	● ● ●	●	● ●	-
Single cycle kinetics	Yes	-	Yes	Plus Package	-	-
Concentration measurement	● ●	● ●	-	● ● Plus Package	●	● ● ●
Calibration-free concentration analysis	Yes	-	-	-	-	-
LMW interaction analysis	● ● ●	● ● ●	-	●	●	-
Thermodynamic characterization	● ● ●	●	-	● Plus Package	●	-
Sample recovery for MS	● ●	-	-	-	● ● ●	-
Performance and technical specifications						
Detection spots/sensor surface	4	20	400	2	4	4
Throughput	● ●	● ● ●	● ● ●	●	● ● ●	●
Unattended run capacity	● ●	● ● ●	●	●	● ●	● ●
Automated data evaluation	● ●	● ● ●	● ●	● ●	●	● ● ●
User guidance	● ●	● ●	● ●	● ● ●	● ●	● ●
Cooled sample storage	● ● ●	● ● ●	-	-	● ●	● ●
Analysis temperature (°C)	4-45 (in-line degasser)	4-40 (in-line degasser)	20-37 (in-line degasser)	25 4-40 with Plus Package, (including in-line degasser)	4-40	25
GxP compliance support	GxP Package	GxP Package	-	-	GxP Package	Yes

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





Mix

Flow/Stop Flow

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<i>Command</i>	<i>Comments</i>	<i>Injected volume</i>	<i>Sample consumption</i>	<i>Max flow rate</i>
INJECT	General purpose inject.	5-325 µl	Injected vol +30 µl	100 µl/min
KINJECT	Low dispersion inject for kinetics and for samples sensitive to dilution with running buffer. Controlled dissociation time after sample injection.	10-250 µl	Injected vol +40 µl	100 µl/min
QUICKINJECT	Low sample consumption. Less preparation time than INJECT but more sample dispersion.	5-325 µl	Injected vol +10 µl	100 µl/min
COINJECT	Injection of second sample immediately following the first (e.g. for dissociation in the presence of ligand). Variable contact times for first sample only.	Sample 1: 10-100 µl Sample 2: 15-250 µl	Sample 1: injected vol +40 µl Sample 2: injected vol +40 µl	100 µl/min
BIGINJECT	Large volume inject. Direct inject up to 325 µl, sample loop between 326 and 425 µl, direct between 426 and 750 µl.	326-750 µl	Injected vol +52 µl	50 µl/min*
MANUAL INJECT	Manually controlled inject. See Section 4.3.9.	1-100 µl / loading	Loaded vol +30 µl	100 µl/min

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- **Cancer and Proteomics Booklet**
- **Neuroscience Booklet**
- **Immunology Booklet**
- **Biacore™ biosensor assays for quantitation of influenza virus and HCP**
- **Hit validation using Biacore™ and MicroCal™ analysis**
- **Unraveling the mechanisms of RNA-binding protein functions using SPR-based kinetic analysis**
- **Initial studies towards automated screening for enzyme inhibitors by SPR-MS**
- **Rapid development of a GMP-compliant assay for the determination of antibody concentration**
- **Validation of a concentration assay using Biacore C**
- **Optimizing SPR-MS for isolation and characterization of low level impurities in complex, protein pharmaceutical samples**
- **Development of immunotherapeutics and immunization regimes**
- **Transition state thermodynamic analysis using Biacore T100**
- **Selection of protein therapeutic candidates using Biacore T100**
- **A multi-protein panel approach for selectivity-based compound screening using Biacore A100**
- **Early kinetic screening of hybridomas for confident antibody selection using Biacore A100**
- **The integration of a protein interaction array into a fully automated selection and screening process**
- **Resolving a bottleneck in screening and characterization of recombinant antibody fragments using Biacore 4000-AA**
- **Analyte recovery in Biacore 3000: optimized functions for SPR-MS applications**
- **Label-free interaction analysis in real-time using surface plasmon resonance**
- **Developing a natural partnership for rapid, sensitive identification of binding partners**

You are here: [Biacore Life Sciences](#)

Interaction Analysis

- [» Introduction](#)
- [» LMV Drug discovery](#)
- [» Biotherapeutics](#)
- [» Research](#)

About Biacore Systems

- [» Introduction](#)
- [» Product Range](#)
- [» Technology](#)
- [» History](#)

Events & Press Room

- [» Events](#)
- [» Webinar Recordings](#)
- [» Press Room](#)

Support & Training

- [» Application Support](#)
- [» Training](#)
- [» Product Support & Service](#)

Web Services

- [» On-line Purchasing](#)
- [» Contact & Ordering](#)
- [» Request Information](#)
- [» My Account](#)

Website Destinations

- [» High Content Analysis](#)

Biacore

label-free interaction analysis –
from research through drug discovery
and development to QC

» Events


- Jul. 13-14, 2010
[BITC Workshop](#)
- Aug. 1-5, 2010
[Protein Society 2010](#)
- Aug. 8 - Jul. 14, 2010
[GRC Medicinal Chemistry](#)
- Sep. 5-9, 2010
[EFMC-ISMIC 2010](#)
- Sep. 19-23, 2010
[HUPO 2010](#)

» Login

Get access to valuable additional material, such as detailed product information, support tools and more.

- [» Sign-up](#)
- [» Benefits of signing-up](#)
- [» Login](#)

Human Fab Capture Kit

Biacore™ 4000

Visit our Training Portal!
A comprehensive range
of learning opportunities

[»MORE INFORMATION](#)

Developments in Protein Interaction
Analysis - Barcelona, October 2010
Early registration discount!

[»MORE INFORMATION](#)