


## 1. Chemical kinetics

### 1.a. Measurement of sucrose inversion in acidic solution



The decomposition of sucrose into glucose and fructose in water (so-called inversion of sucrose) is a reaction that can be accelerated by acidic catalysis:



In excess of water, it is a pseudo-first order reaction and the rate equation is given by the expression:

$$v = -\frac{dc}{dt} = k \cdot c \quad (1.2.)$$

where  $c$  is the sucrose concentration,  $t$  means time, and  $k$  is the velocity constant of the acid catalysed reaction. With respect to the catalyst concentration, the following applies:

$$k = k_0 + k_1[H^+] \quad (1.3.)$$

where  $k_0$  is the rate constant in pure solvent ( $H_2O$ ) that is practically negligible,  $k_1$  is the rate constant of the reaction that is homogeneously catalysed by hydrogen ions of concentration  $[H^+]$ .

The first-order rate constant is given by the relationship:

$$k = \frac{1}{t} \ln \frac{c_0}{c} \quad (1.4.)$$

where  $c$  is sucrose concentration at time  $t$  from the reaction beginning (ie from acidification of the solution),  $c_0$  is the initial sucrose concentration.

The sucrose solution turns the plane of polarized light clockwise while the equimolar mixture of the fructose and glucose turns the light counter-clockwise after complete reaction. We can detect a shift from initial angle  $\alpha_0$  of polarized light to final angle  $\alpha_\infty$ .

The fraction  $c_0/c$  in relation (1.4.) can be replaced by a fraction, which depends on the angle rotation of the plane of light  $\alpha_t$  at time  $t$ :

We obtain relationship:

$$k = \frac{1}{t} \ln \frac{\alpha_0 - \alpha_\infty}{\alpha_t - \alpha_\infty} \quad (1.5.)$$

where  $\alpha_\infty$  can be determined either after a few days when the reaction is completed, or by computing using empirical equation:

$$\alpha_\infty = -[0.20695 + 0.00004 \cdot c_0 - 0.0032(\tau - 20)] \cdot c_0 \cdot l \quad (1.6.)$$

where  $c_0$  is given in unit  $g/dm^3$ ,  $\tau$  is the solution temperature in  $^\circ C$  and  $l$  is the polarimetric tube length in meters.

The rearranging of eqn (1.5.) gives:

$$t = \frac{1}{k} \ln(\alpha_0 - \alpha_\infty) - \frac{1}{k} \ln(\alpha_t - \alpha_\infty) \quad (1.7.)$$

which is a linear dependence  $t$  on  $\ln(\alpha_0 - \alpha_\infty)$  that allows to determine the rate constant  $k$  from the slope of the regression line.



**TASK:** Measure the turning of the polarized light during the inversion reaction of sucrose and evaluate the rate constant at two different concentrations of the catalyst.



**LABORATORY AIDS AND CHEMICALS:** polarimeter, 2 polarimetric cells, weighting boat, stopwatch, 2 volume pipettes ( $20\text{ cm}^3$ ), volumetric flask ( $100\text{ cm}^3$ ), 3 Erlenmeyer flasks ( $100\text{ cm}^3$ ), sucrose,  $4\text{ M HCl}$  and  $2\text{ M HCl}$ .




**INSTRUCTIONS:** Switch polarimeter on and get to know with instrument operating manual.

Prepare  $100\text{ ml}$  of stock solution by dissolving  $20\text{ g}$  (accurate to  $0.05\text{ g}$ ) of sucrose in water (use  $100\text{ cm}^3$  volumetric flask).

Take in Erlenmeyer flask and mix  $20\text{ cm}^3$  of sucrose stock solution with  $20\text{ cm}^3$  of water and mix thoroughly. Fill the polarimetric cuvette with the solution. Insert the cuvette into the polarimeter (see scheme in **FIGURE 1**) and measure the starting angle of the plane of polarized light for the sucrose solution. Empty the polarimetric cuvette and rinse it with distilled water.

Prepare the reaction solution in Erlenmeyer flask by mixing  $20\text{ cm}^3$  of stock solution of sucrose and  $20\text{ cm}^3$  of  $4\text{ M HCl}$ . Mix the solutions quickly and intensively. Start the stopwatch immediately. Rinse out the polarimetric cell with small amount of reaction solution. Finally, pour the reaction solution up to cuvette marking and do measuring the plane of the polarised plane of light in 5 minute intervals within one hour. Measure the angle also at time 80, 100 and 120 minutes. Similarly prepare a reaction solution applying  $2\text{ M HCl}$  and use same procedure as before. The rotation changes of both reaction mixtures can be measured simultaneously at shifted times.



**FIGURE 1:** Scheme of polarimetric measurement. L-sodium lamp, P-polarizer, A-analyser, K-cuvette with sucrose solution.

**REPORT:** weight of sucrose, value  $\alpha_\infty$  and temperature  $\tau$  used for evaluation of it. **Table 1a** and **1b:** for both experiments at different  $[\text{H}^+]$  concentration: time  $t$ , angle  $\alpha_t$ ,  $\ln(\alpha_t - \alpha_\infty)$ . **Common graph 1:** for both  $[\text{H}^+]$  concentrations: time dependency of angle  $\alpha_t$ . **Common graph 2:** for both  $[\text{H}^+]$  concentrations: dependency  $t$  on  $\ln(\alpha_t - \alpha_\infty)$ . **Next:** Evaluation of rate constants  $k$  at variable  $[\text{H}^+]$  and their confidence intervals. Evaluation of rate constants  $k_0$  and  $k_1$  by solution of two equations with two unknowns (use eqn (1.3)).