

Sequence	Order	Sample	[Product] (mM)	[Enzyme] (mg/mL)
	1	Buffer	0	-
	2	Cal-1	1,5	-
	3	Cal-2	1,2	-
	4	Cal-3	0,9	-
	5	Cal-4	0,6	-
	6	Cal-5	0,3	-
	7	Buffer	0	-
	8	Enzyme 1	To be calculated	0,12
	9	Buffer	0	-
	10	Enzyme 2	To be calculated	0,2
	11	Buffer	0	-

Detected peaks

	Buffer	Cal-1	Cal-2	Cal-3
	1,6012	0,4512	0,8517	1,1097
	1,6040	0,4287	0,8151	1,0760
	1,6063	0,4411	0,7821	1,0599
	1,6011	0,4493	0,8157	1,0732
	1,5987	0,4566	0,7959	1,0846
	1,6203	0,4503	0,7834	1,1003
	1,6059	0,4472	0,7826	1,1002
	1,6504	0,4413	0,7839	1,0830
	1,6096	0,4356	0,7659	1,0869
	1,6154	0,4303	0,7748	1,0754
Average				
Product (mM)				
Activity ($\mu\text{mol}\cdot\text{s}^{-1}\cdot\text{mg}^{-1}$)				

Parameters

Reaction time

4 min

Visualized data

Cal-4	Cal-5	Buffer	Enzyme 1	Buffer	Enzyme 2	Buffer
1,3497	1,5202	1,6581	1,0675	1,5084	0,7857	1,4512
1,2778	1,4962	1,6310	1,0096	1,5539	0,7707	1,5199
1,3116	1,4665	1,6435	0,9797	1,5499	0,7426	1,5520
1,2890	1,4579	1,6354	0,9770	1,5935	0,7532	1,5669
1,2337	1,4831	1,5797	0,9738	1,6050	0,8608	1,5978
1,3331	1,4690	1,6503	0,9648	1,6444	0,7062	1,6033
1,3396	1,4789	1,5723	0,9662	1,6214	0,7087	1,5885
1,3509	1,4785	1,5921	0,9642	1,6222	0,7171	1,6017
1,3090	1,4653	1,6599	0,9574	1,6501	0,7142	1,6049
1,3268	1,4587	1,6448	0,9364	1,6198	0,7086	1,6223

--

Procedure

1. Calculate the average

Note: Calculate the average for each sample.

2. Construct the calibration curve (dependence of signal on concentration) with a linear regression curve.

Note: Include also the blank concentration.

3. Determine [P] de Enzyme 2.

4. Based on this info calculate the specific activities at each concentration using the following formula:

$$a = \frac{[P]}{[E] \cdot t}$$

a = specific activity
[P] = product concentration
t = reaction time
[E] = enzyme concentration

Note: the specific activity is expressed in units of $\mu\text{mol} \cdot \text{s}^{-1}$ which units do you use?

5. Compare both enzymes

average of each sample.

average of all buffers as if it

calibration curve (the
relationship on [P] - product
either linear or quadratic

the buffers with zero product

developed by Enzyme 1 &

information, calculate their
according to the following

concentration

concentration

activity should be expressed
[$\mu\text{mol}\cdot\text{mg}^{-1}\cdot\text{s}^{-1}$], so double-check
input into the formula.

enzymes.