
MASARYKOVA UNIVERZITA

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<https://www.sartorius.com/en/product/product-detail/16757/>

BioSampler

8-hour Sampling of Bioaerosols into Liquid

- Patented[®] design overcomes problems associated with impinger sampling of bioaerosols
- Allows use of non-evaporating collection liquids for longer sampling times
- Maintains constant sampling efficiency up to 8 hours
- Patented[®] swirling liquid collection method
 - Significantly reduces particle bounce and re-aerosolization
 - Preserves microorganism integrity and viability
- Collection liquid easily transferred to agar plate for culturing
- Samples can be analyzed by a variety of methods
- Reusable — can be autoclaved

The SKC BioSampler[®] is a glass collection device that externally resembles an All-glass Impinger (AGI-309). Internally, BioSampler contains design features that overcome some of the sampling problems associated with using impingers for bioaerosol collection.

- Inlet limits collection of airborne particles to those that would pass through the human nose.
- Three tangential nozzles reduce particle bounce off inner wall.
- Airflow through the nozzles causes the collection liquid to swirl upward on the

The BioSampler Advantage

- **Constant sampling efficiency over longer sampling times**
 BioSampler: Thicker, non-evaporating liquids such as ViaTrap[®] mineral oil can be used to maintain constant sampling efficiency over an 8-hour workshift. Longer sampling times increase sample volumes for detecting organisms at lower concentration levels.
 Standard Impingers: Typical sampling times with standard impinger liquids are only 1 to 1.5 hours.

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