

Solvent Suppression using TopSpin 3.x

updated: 15 Jan 2013 (cgf)

Introduction

There are number of solvent suppression schemes used in NMR, and the technical details can become extensive. Here is a brief list [the two- (or more) letter code identifier in par and pp names]:


- presaturation** – often gives cleanest suppression, but removes exchangeable protons [**pr** or **ps**]
- noesy1d preset** – can provide cleaner baselines; used in many metabolomics studies [**pn**]
- purge preset** – a relatively new method of presaturation, often superior to standard preset [**purge**; **pp** is a purge pulse that is different]
- wet** – a prominent method used with mixed solvents (multiple-peak suppression); also reduces intensity of exchangeable protons [**wt** or **wet**]
- watergate** – various types exist, the basic requiring optimization but is high-quality [**wg**]
- 3-9-19** – the most common wg variant; simple to optimize with no effect on exchangeable protons [**wg** or **19**]
- w5** – a relatively new wg variant similar to 3-9-19 [**w5**]
- excitation sculpting** – a high-quality method utilizing selective rf and gradients [**es**]

A significant issue with solvent suppression on Bruker spectrometers is which variants are available in the experiment you really need. **purge** may run great in a ^1H 1d, but does not exist (currently) in tocsy or noesy 2Ds. All the following describe only how to obtain 1D ^1H spectra. See cgf for more information about extension into 2D or other experiments. *Another issue is quantitative behavior. Presat, wet and watergate are good quantitatively (except for exchangeable protons with presat). Do not use noesy1d quantitatively. Until the shaped-rf is updated, do not use exct sculpting either.*

A. Initial Setup:

1. Always start by acquiring a one-scan ^1H spectrum.

It is not required, but best to run solvent suppression experiments on-resonance to the solvent peak. If you believe this is not optimal, find cgfry for further discussion.

2. Put the solvent peak on-resonance by clicking , left-click with the cursor in the middle of the solvent peak, then choose O1. Retake the one-scan ^1H spectrum to verify that the peak is in the center of the spectrum.

[Another, possibly simpler, method for going on resonance is to enter the shift of the solvent peak in the parameter **o1p**. This will not be as accurate, since it won't account for solution/temperature/etc variations in chemical shift.]

3. For many methods — watergate, wet, excitation sculpting — the above is sufficient. But for preset methods, it is best to optimize the **o1** (**o1p**) value by using the **gs** utility:

- a) **gs**↵
- b) click **offset**
- c) change **o1** by click on sliding scale above/below; wait 1 to 2 scans to see effect
- d) SAVE on best value; for most sequences here, **o1** (not **o2**, definitely not **o3**)

4. Copy parameter into a new expno using: **iexpno**

B. Presaturation:

1. In the new expno, run **ased** and change the first parameter PULPROG to **zgpr**.
2. The critical new parameter is **PLW9**, which will perform a low-power cw pulse on-resonance. You can raise the power of this parameter to decrease the intensity of the residual signal, but setting it too high will damage the probe!

$$\begin{aligned} \text{PLW9} &\leq 3 \text{ mWatts } (\leq 0.003 \text{ in the 1}^{\text{st}} \text{ box on the ased screen}) \\ &\geq 25 \text{ -dBW } (\text{value in 2}^{\text{nd}} \text{ box on the ased screen should } \geq 25) \end{aligned}$$

3. Run **rga** prior to doing **zg**.
4. Take new data with **ns = 4×i** and **ds = 2×j** (i,j = 1,2,4,8,...).
5. **zgcpr** and **zgcpgr** are two variants using composite pulse. So, cgf's experience does not recommend these (but **zgcpgr** is the better of the two).

C. noesygpr1d Presaturation:

1. Do the same steps above, but read in the Bruker parameter set **watersup**. I.e., do
rpar WATERSUP all
2. Change **o1** or **o1p** to match the value found in A.2.
3. Run **rga** prior to doing **zg**.
4. Take data with **ns = 8×i** and **ds = 4×j** (i,j = 1,2,4,8,...).
5. PLW9 again might be smaller than optimal. Same conditions apply as in B.2.

D. purge Presaturation:

1. Do the same steps above, but read in the UWchem parameter set **probe_1Hpurge**. I.e., do
rpar probe_1Hpurge all
2. Change **o1** or **o1p** to match the value found in A.2.
3. Run **rga** prior to doing **zg**.
4. Take data with **ns = 1×i** and **ds = 2×j** (i,j = 1,2,4,8,...).

E. watergate 3-9-19 suppression:

1. Do the same steps above, but read in the Bruker parameter set **P3919GP**. I.e., do
rpar P3919GP all
2. Change **o1** or **o1p** to match the value found in A.2.
3. Change **d19 = 180µs**. The spectrum will have “nulls” in peak intensities (excitation profile) at $1/(2 \times d19)$, so this is a critical parameter. 180µs works well at 500 MHz (nulls at ± 5.5 ppm from spectrum center).
4. Run **rga** prior to doing **zg**.

5. Take data with **ns = 8×i** and **ds = 4×j** (i,j = 1,2,4,8,...).

F. watergate w5 suppression:

1. Do the same steps above, but read in the Bruker parameter set **P3919GP** . I.e., do
rpar P3919GP all
2. Run **ased** and change the first parameter PULPROG to zgpgw5.
2. Change **o1** or **o1p** to match the value found in A.2.
3. Change **d19 = 180µs** [see comments in section E.3 above]
4. Run **rga** prior to doing **zg** .
5. Take data with **ns = 8×i** and **ds = 4×j** (i,j = 1,2,4,8,...).

G. excitation sculpting:

1. Do the same steps above, but read in the Bruker parameter set **P3919GP** . I.e., do
rpar ZGESGP all
2. Change **o1** or **o1p** to match the value found in A.2.
3. Have to figure out how to get away from the rectangular selection pulses (these will produce “weird” peak intensities rippling throughout the spectrum)!!
4. Run **rga** prior to doing **zg** .
5. Take data with **ns = 8×i** and **ds = 4×j** (i,j = 1,2,4,8,...).

H. Multiple-peak solvent suppression:

1. Take a 1-scan proton spectrum as described in section A.
2. Integrate the solvent peaks to be suppressed in the proton spectrum, and **Save As... reg** : see section B of the notes at:
http://www.chem.wisc.edu/~cic/nmr/Guides/Ba3vug/HW637/HW5_Av400-sel1DbyTopSpin.pdf
3. click on: **CREATE DATASETS**
and choose either preset or wet toward the bottom of the list.