

## Topic P10: Fundamentals of clinical mycology

To study: Fungi

From spring term: Microscopy, culture, antibiotic susceptibility, precipitation

### Table for major results of Task 1 to Task 3 (to be filled step by step):

Strain		K	L	M	N
Gram stain – Task 1					
Culture (blood agar) Task 2a	Size				
	Colour				
	Shape				
	Profile				
	Haemolysis				
	Surface				
	Odour				
Task 2b: Growth on Sabouraud + chloramphenicol (growth +/-)					
Task 2c: Growth on a chromogenic medium					
<b>FINAL CONCLUSION (according to Task 2c/Task 3)</b>					

### Task 1: Microscopy of bacteria and yeasts

Stain the presented cultures of microorganisms by Gram. Use immersion microscopy (immersion objective magnifying 100×). Write your results in the table. Note the size difference between yeasts and bacteria.

### Task 2: Culture of bacteria and yeasts

#### a) Culture on blood agar

Describe colonies of the presented strains on blood agar and fill in the main table. Do not forget to describe the odour. Note that the colonies of yeasts (according to Task 1) are similar to colonies of some bacteria (especially to the colonies of G+ cocci, according to the morphology of colonies most probably staphylococci).

#### b) Culture on Sabouraud agar with chloramphenicol

Evaluate growth of the presented strains on selective agar for yeasts and moulds (Sabouraud agar with chloramphenicol). Sabouraud agar itself is not selective but it is made selective for yeasts and moulds using broad-spectrum antibiotics (chloramphenicol).

#### c) Culture on a chromogenic medium

Chromogenic media for yeasts enable differentiation of the most important species of the *Candida* genus. With the help of control strains, try to determine the *Candida* species using the chromogenic medium. If the colonies are white (no colour substance present) the strain cannot be identified using this chromogenic medium.

### Task 3: Biochemical identification of yeasts using Auxacolor set

Read the results of Auxacolor for the strain that could not be identified using the chromogenic medium. In this case, we do not calculate the code but we compare the results with the table. (You might obtain more than one taxon. If so, remember that *C. albicans* would be green in 2c) and *Rhodotorula* would be pigmented red in 2b.)

C Neg*	Glu**	Mal**	Sac**	Gal**	Lac**	Raf**	Ino**
Cel**	Tre**	Ado**	MeI**	Xyl**	Ara**	Act***	Pox#
Result of identification:							

\*normal blue

\*\*\*yellow positive, colourless negative

\*\*yellow positive, blue negative

#brown positive, colourless negative

**Task 4: Assessment of antimicrobial drugs susceptibility**

For the treatment of fungal infections, it is not possible to use common antibiotics. We have to use special drugs – antimycotics. These, on the other hand, are not effective against bacteria.

**a) Assessment of susceptibility to antimycotics**

Perform in vitro susceptibility testing of the presented strains to antimycotics. Use the table only for two strains that are susceptible to at least something, the other two strains should be written beneath the table. Complete the table with the full name of the antimycotics according to the card and for all the tested strains assess the susceptibility or resistance. For Amphotericin B, the reference zone is 10 mm. For other antimycotics, it is 20 mm, but the inside of the zone does not have to be perfectly clear.

Strain →				
Antimycotics (full name)	Zone Ø (mm)	Interpretation	Zone Ø (mm)	Interpretation

Strains \_\_\_ and \_\_\_ (i. e. strains of bacteria–fungi\*) are resistant to all the tested drugs. *\*delete as appropriate*

**b) Assessment of susceptibility to antibiotics**




Perform in vitro susceptibility testing of the presented strains to antibiotics. Use the table only for two strains that are susceptible to at least something, the other two strains should be written beneath the table. Complete the table with the abbreviations of the antibiotics according to the card and for all the tested strains, measure the susceptibility zones. On your card there are limit zones – according to these, interpret the zones as susceptible, (S) resistant (R) and dubious (D).

Strain →				
Antibiotic (full name)	Zone Ø (mm)	Interpretation	Zone Ø (mm)	Interpretation

Strains \_\_\_ and \_\_\_ (i.e. strains of bacteria–fungi\*) are resistant to all the tested drugs. *\*delete as appropriate*

**Task 5: Microscopy of moulds**

Moulds are usually microscoped differently than yeasts. Gram staining is used rarely. **Mostly we use wet mount, objectives magnifying 10× to 40× (→ NO IMMERSION OIL!).** Draw and describe the three presented species of moulds (if they are more, choose three of them). Do not forget to describe hyphae, macro- and microconidia and other observed objects. Use lines or numbers to connect the description with the objects.


Genus	Genus	Genus
Species	Species	Species
		
Description:	Description:	Description:
_____	_____	_____
_____	_____	_____

**Task 6: Culture of moulds**

Moulds usually require longer time to grow (weeks). That is why we usually do not use Petri dishes for culture but only test tubes (to avoid drying and contamination). Draw culture results of the presented moulds. (If they are more than 3, choose 3 of them).

**Task 7: Indirect diagnostics of aspergillosis**

Evaluate the results of precipitation in gel in order to diagnose antibodies in aspergillosis. Draw the result and mark the positive and negative cases.



drawing  
description

**Task 8: Sampling for mycoses**

Watch the video "Sampling of cutaneous mycoses". Write down basic rules for sampling in superficial mycoses.

**Attention! For the next practical, one volunteer for each group is needed. The volunteer should not clean his/her teeth in the morning. It will be possible to clean the teeth immediately after the practical, and there will be no permanent effects on your teeth.**