PZ13 Clinical microbiology IV – examination of wound and bloodstream infections

To study: Your own protocols (especially Special bacteriology)

Wound infections

Task 1: Specimens in wound infections

Try to fill in the following table:

Type of wound	Superficial	Deep wound with	Deep wound with not	Wound with pus,		
	wound	amount of pus sufficient	sufficient amount of	possibly containing		
		for being sent as a liquid	pus	anaerobic bacteria		
Sampling method						
When a specimen from a wound is send to the laboratory, it is very important to fill in the request form,						
especially to write 1) and 2)						

Task 2: Imprint method for superficial wound examination (moulage method)

a) Imprint method – performing

A sterile filtration paper on is placed on a superficial wound. We let it for 10 seconds here, then using tweezers, we transport it carefully to a Petri dish with nutrient agar. After that, the filtration paper is sent together with the agar plate to the laboratory. In the laboratory the filtration paper is placed to two or three more media: agar with 10 % NaCl, chromogenic URI medium etc. After that, all media are cultivated overnight. Dental students do not perform this part practically.

b) Imprint method – reading of results

Try to read the result of imprint method on URIchrom chromogenic medium using recounting scheme on your table and with the help of the key of colours of individual bacteria on the chromogenic medium. Attention! You have real results from real patients. Your result is not supposed to be the same as the result of your neighbour with another agar plate. Even the number of strains may be different.

The cultivation result of my imprint contained:

Likely species of bacterium	Quantity (approx. number of colonies per 25 cm ²)
1.	
(2.)	
(3.)	

Clue for preliminary diagnostics: Staphylococci – white on URI, growing also on NACL, white colonies on blood agar; Haemolytic streptococci – haemolytic colonies on blood agar, not growing on NACL, on URI not growing or (S. agalactiae) pale blue. Enterococci have grayish colonies on URI and small, but clearly blue colonies on URI. Enterobacteriaceae and G- non-fermenters – growing on Endo agar. Escherichia is pink on URI, Klebsiella is blue on URI, Proteus is yellow on URI, Pseudomonas is white or slightly green (because of its own pigmentation) on URI. All this is only preliminary, the algorithms from previous practicals are valid!

Task 3: Deeper wound swab result

In the case of a wound swab, there is no "common flora". That is the main difference between wound swab and e. g. swabs from respiratory ways: it is not necessary to search for a pathogen among the normal flora. On the other hand, we mostly use more culture media to detect all possible pathogens, even if they would be in a mix of them. Besides blood agar and Endo agar we usually use also blood agar with 10 % NaCl and blood agar with amikacin in order to search for streptococci and enterococci (but none of these media is used in our task). In other situations there is one pathogen only, and even in small amounts, so we have to multiply it in a liquid medium (broth). Also this medium is not present in our task. Fill in the form again.

Name	Dental Medicine	Date 9. 12.	2013	Page 1/4
Name	Deniai Medicine	Date 9. 12.	2013	1/4 Tage 1/4

Č. pojištěnce *1983 Dg: Suppurating wound of Val	provedl	9	1 - 1
Pacient Lucy Yellow C. pojištěnce *1983 Variabilní symbol Odeslán ad: Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: Poznámka: Dg: Suppurating wound of planta pedis Kód náhrady Požadováno: Kód náhrady Poznámka: To Dr. Micube Tephole 10 Dne: Dne: Dne: Pod Suppurating wound of planta pedis Val Val Val Val Val Val Val Va	ALTERNATION OF THE PARTY OF THE	Poř. č.	7
Pacient Lucy Yellow Ĉ. pojištěnca *1983 Variabilní symbol Odeslán ad: Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; The pus appeared after two days Poznámka: Poznámka: Dg: Suppurating wound of planta pedis Kód náhrady Požadováno: A Poznámka: Dne: Dne: Dne: Dne: Doda Do	រាវា គ្រ		
C. pojištěnce *1983 Variabilní symbol Odeslán ad: Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: Dg: Suppurating wound of planta pedis Kód náhrady Z Ród náhrady Z Poznámka: Da Odeslán ad: I Kód náhrady Z Poznámka: Da Odeslán ad: I Kód náhrady Z Da Odeslán ad: I Kód náhrady Z Da Odeslán ad: I Da Odeslán ad	ČP		1. 17
Variabilni symbol Odeslán ad: Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: Poznámka: T2 Dr. Microbe Termole 10	Odbornost	4	of Fig
Odeslán ad: Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: To Dr. Micube Terrole 10 123 general ractif oner Dne: 11	/ar. symbol		
Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: T2 Dr. Micube Termole 10 Dne: 11	Datum	Kód	Poč
Požadováno: Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: 72 Dr. Micube Terrible 10 Dne: 11			
Wound with pus on planta pedis, caused by stepping on a tin in a pond; the pus appeared after two days Poznámka: 72 Dr. Microbe Termole 10 Dne: Dne: 11			
by stepping on a tin in a pond; the pus appeared after two days Poznámka: To Dr. Miculae Terrible 123 general rescrit oner Done:			
by stepping on a tin in a pond; he pus appeared after two days Poznámka: 72 Dr. Miculos Terrible 123 general rescrit oner Dne: 11		1	
Poznámka: Tamber Terrible 123 general kractit oner The pus appeared after two days 7 7 8 9 10 Dne: 11			
Poznámka: 7	1 1 1		
72 Dr. Micurbe Teyrole 10 123 general reactif oner Dne:			18.5
72 Dr. Micurbe Teyrole 10 123 general reactif oner Dne:			
123 general ractif oner Dne:			1
123 general ractif oner Dne:			100
		بننث	1 2
456 Grampositive 8, Brno			1.
azítko a podpis léko			
VZP-08x1999 razitko a podpis	(1 (C) (4)		

Patient:Lucy	84 Dg∴wound of planta pedis						
Specimen: wound swab* Ordered by: Dr. Microbe Terrible							
*note: pyogene wound on planta pedis, swimming in a pond							
Growth on blood a. (incl. sme	ll) Endo agar:	MH agar:		Oxidase:	Conclusion:	Interpretation	
Antibiotic susceptibility test				al conclusion treatment:	and recommen	dment	

Name ______ Dental Medicine Date 9. 12. 2013 Page 2/4

Bloodstream infections

Task 4:	Blood	cultures –	processing
I ask T.	Diou	cuitui cs —	DI OCCSSIIIZ

Describe the use of three types of blood culture vessels.
Fill in which data should not be missing on the order form in the case of blood culture (only "materia"
type/examination type" field)
Familia.
Explain: Why is absolute sterility in blood culture samples more necessary than in any other blood specimens (e. g. those
sent for biochemical examination)?
II
How many blood cultures should be taken and why?
Fill in the missing fields in the description of blood culture processing and examination according to the vide clip and the teacher's explanation.
A blood culture vessel arrives in the laboratory. Here it is put into a
The positive result is demonstrated by and
When the cultivation is positive, a smear is prepared and the content of the vessel is
onto the blood and Endo agar. Also, a preliminary test is performed directly
from the specimen; as the inoculum is not standardized here, its results are only
Task 5: Blood cultures – microscopy of a positive specimen
The cultivator for blood cultures revealed a positive result. For preliminary treatment, a Gram stained smear
performed from the contain of the vessel. Observe the result and write it. Attention! The slides have origin it real blood cultures of different patients. Therefore your result may be simply different from that of you
neighbour with a different slide.
Blood culture contained gram-positive – gram-negative* cocci – bacilli* arranged in** * delete as appropriate **only for cocci (pairs, chains, clusters) or G+ bacilli in palisades
Task 6: Blood cultures – cultivation result Observe cultivation result of a positive blood cultures inoculated on solid media. Suggest more methods for
detailed diagnostics of bacteria. Try to assess preliminary antibiotic susceptibility. Also here you are no
supposed to have the same results as your neighbour.
Name of medium Crowth V/N appearance
Growth Y/N, appearance of colonies

Name _____ Dental Medicine Date 9. 12. 2013 Page 3/4

More tests of more detailed determination:							
Preliminary name of the microbe: Susceptibility testing Name of the set of antibiotics:							
Antibiotic	Reference size	Measured size	Susceptible - resistant	Antibiotic	Reference size	Measured size	Susceptible - resistant
1.			S-R	4.			S-R
2.			S-R	5.			S-R
3.			S-R	6.			S-R

Task 7: Blood cultures – interpretation

Look at interpretation for results of two different patients.

John White, *1942, elevated temperature and	Joe Black, *1945, elevated temperature and			
inflammatory markers, three blood culture specimens	inflammatory markers, three blood culture specimens			
sent to the laboratory	sent to the laboratory			
I Central venous catether. Time to detection 10 hours,	I Central venous catether. Time to detection 8 hours,			
finding: Staphylococcus hominis, susceptible to	finding: Staphylococcus epidermidis, susceptible to			
oxacilin, tetracycline, vankomycin, resistant to	oxacilin, resistant to tetracycline, vankomycin,			
erythromycin, klindamycin, co-trimoxazole.	erythromycin, klindamycin, co-trimoxazole.			
II Peripherial catather. Time to detection 13 hours,	II Peripherial catather. Time to detection 26 hours,			
finding: Staphylococcus hominis, susceptible to	finding: Staphylococcus hominis, susceptible to			
oxacilin, tetracycline, vankomycin, resistant to	oxacilin, tetracycline, vankomycin, erythromycin,			
erythromycin, clindamycin, co-trimoxazole.	clindamycin, co-trimoxazole, no resistance observed			
III Venepunction. Time to detection 13.5 hours,	III Venepunction. Time to detection 38 hours, finding:			
finding: Staphylococcus hominis, susceptible to	Staphylococcus epidermidis, susceptible to oxacilin,			
oxacilin, tetracycline, vankomycin, resistant to	co-trimoxazole, vankomycin, resistant to tetracycline,			
erythromycin, clindamycin, co-trimoxazole.	erythromycin, clindamycin.			
Likely interpretation:	Likely interpretation:			
Probably bacteriaemia	Probably pseudobacteriaemia			
Provavsy vacieriaemia	Provavvy pseudovacieriaemia			

Name ______ Dental Medicine Date 9. 12. 2013 Page 4/4