

Lab class no. 2 – Evaluation of efficacy of disinfectants, evaluation of colonies, Gram staining, isolation inoculation

(in pairs)

Aims:

- *Evaluation of efficacy of disinfectants.*
- *Evaluation of colonies*
- *Gram staining*
- *Microscopy – permanent preparations*

Material

Plastic bacteriological loop

Set of solutions for Gram staining

Glass slides

Ethanol burner

Microscope

Prepared dish Nutrient broth agar and Malt agar for isolation inoculation

Work procedure

1. Evaluation of efficacy of disinfectants and evaluation of colonies of microorganisms *M. luteus*, *E. coli* and *S. cerevisiae*:

- count the colonies in the respective quadrants and calculate % of efficacy of all disinfectants for all three microorganisms - put the results of % of efficacy to graphs into your protocols
- further evaluate colonies of all three microorganisms with respect to these parameters:
 - a. **size** - colonies are dotted (if smaller than 1 mm) or their size is expressed as diameter in mm
 - b. **shape** - round, irregular shape, ...
 - c. **edges** - straight, wavy, protuberant, ...
 - d. **profile** - flat, slightly bulging, convex, raised edges, ...
 - e. **surface** – shiny (S-phase), smooth, opaque, rough (R-phase), wrinkled, ...
 - f. **transparency** - transparent, translucent, opaque

- g. **color** - colorless, white/greyish, production of pigments (*M. luteus* – yellow, *M. roseus*, *Serratia rubidea* – rose) or change in color of indicator present in medium (lactose fermenting and non-lactose fermenting enterobacterias on selective diagnostic media with lactose)
- h. **change in surroundings** - change in color (production of pigments soluble in water, change in media pH – change in color of indicator), on blood agar we can see hemolysis of erythrocytes (partial hemolysis - α , full - β , double - α' , greening - viridation)

2. Gram staining of *M. luteus*, *E. coli* and *S. cerevisiae*:

- Transfer small amount of chosen colony with bacteriological loop to a drop of physiological solution on glass slide and let dry
- this preparation then fix using heat - three times pass the glass slide through a flame; reversed side facing the flame!! (mild denaturation and killing of bacteria)
- now add to the preparation few drops of crystal violet for 20 seconds
- without discarding, add few drops of Lugol solution
- discard both solution and add enough Lugol solution and let it act for 20 - 30 seconds
- using a pipet add acetone-ethanol mixture to glass slide and let it flow over the preparation, until it becomes discolored (mixture flowing from preparation is also discolored)
- preparation then wash with distilled water
- then add few drops of safranin solution and let it dye the preparation for 30 - 60 seconds
- finally, wash the preparation and let dry; microscope in oil immersion (objective 100 \times).

3. Microscopy – permanent preparations

Protocol no. 2 will contain:

- evaluation of colonies of all three microorganisms as described in point no. 1
- procedure of Gram staining, principle, drawing or picture of your own results

- drawing of photo of permanent preparations - describe: cocci, bacilli, G+, G-, aerobic, anaerobic, normal microbiota or not, where we can find them (skin, GIT, ...), disease which they can cause