MASARYKOVA UNIVERZITA

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<http://www.sci.muni.cz/ccm>

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
BioSampler

8-hour Sampling of Bioaerosols into Liquid

- Patented* design overcomes problems associated with impinger sampling of bioaerosols
- Allows use of non-evaporating collection liquids for longer sampling times
- Maintains constant sampling efficiency up to 8 hours
- Patented* swirling liquid collection method
 - Significantly reduces particle bounce and re-aerosolization
 - Preserves microorganism integrity and viability
- Collection liquid easily transferred to agar plate for culturing
- Samples can be analyzed by a variety of methods
- Reusable — can be autoclaved

The SKC BioSampler® is a glass collection device that externally resembles an All-glass Impinger (AGI-30). Internally, BioSampler contains design features that overcome some of the sampling problems associated with using impingers for bioaerosol collection.

- Inlet limits collection of airborne particles to those that would pass through the human nose.
- Three tangential nozzles reduce particle bounce off inner wall.
- Airflow through the nozzles causes the collection liquid to swirl upward on the

The BioSampler Advantage

- **Constant sampling efficiency over longer sampling times**

BioSampler: Thicker, non-evaporating liquids such as Via Trap® mineral oil can be used to maintain constant sampling efficiency over an 8-hour workshift. Longer sampling times increase sample volumes for detecting organisms at lower concentration levels.

Standard Impingers: Typical sampling times with standard impinger liquids are only 1 to 1.5 hours.

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Table 4 : Categories of CFU/m³ (mixed populations of fungi) obtained with the Andersen six-stage sampler in combination with MEA, and with the N6-Andersen one-stage sampler in combination with MEA and DG18.

category	houses (CFU/m ³)	non-industrial indoor environments (CFU/m ³)
very low	< 50	< 25
low	< 200*	< 100
intermediate	< 1 000	< 500
high	< 10 000	< 2 000
very high	> 10 000	> 2 000

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EUR 14988 (Report No. 12: Biological Particles in Indoor Environments, Commission of the European Communities, Report No. 12, Luxembourg, 1994)

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