



Central European Institute of Technology
BRNO | CZECH REPUBLIC

Boris Tichý

DNA re-sequencing


The logo for MUNI (Masaryk University) features the letters "MUNI" in a large, white, sans-serif font, with each letter partially cut off on the right side.



DNA re-sequencing

Comparisson with reference sequence (genome)

- WGS – whole genome sequencing
- WES – whole exome sequencing
- Genes, gene panels
- SNV – single nucleotide variant
- CNV – copy number variant




Nature Reviews Genet. 2010 Oct ;11(10):685-96.

DNA re-sequencing

WGS – whole genome sequencing


- DNA sequencing without selection
- Fragmentation
 - Ultrasound
 - Enzymes
 - Longer fragments (> 300bp)
- PCR-free – GC-bias
- Paired-end sequencing
- Germline - 30x coverage ~ 100Gb



DNA re-sequencing

WES – whole exome sequencing


- Coding sequence selection
- Fragmentation
 - Ultrasound
 - Enzymes
 - Fragments 150-250bp
- PCR
- Paired-end sequencing



DNA re-sequencing

WES – whole exome sequencing


- Coding sequence selection
- Fragmentace
 - Ultrasound
 - Enzymes
 - Fragments 150-250bp
- PCR
- Paired-end sekvenace
- Germline - >50x coverage, >3Gb



DNA re-sequencing

Gene panels


- Coding sequence selection
- Fragmentation
 - Ultrasound
 - Enzymes
 - Fragments 150-250bp
- PCR
- Paired-end sekvenace
- Germline - >50x coverage



DNA re-sequencing

Gene panels


- Coding sequence selection
- **PCR enrichment**
- Single primer PCR
- Inverse PCR + capture
- Capture + probe extension
- Paired-end sequencing
- Germinální - >50x coverage



DNA re-sequencing

Gene panels


- Coding sequence selection
- PCR enrichment
- **Single primer extension**
- Inverzní PCR + capture
- Capture + probe extension
- Paired-end sekvenace
- Germinální - >50x coverage



DNA re-sequencing

Gene panels


- Coding sequence selection
- PCR enrichment
- Single primer PCR
- **Inverse PCR + capture**
- Capture + probe extension
- Paired-end sekvenace
- Germinální - >50x coverage



DNA re-sequencing

Gene panels

- Coding sequence selection
- PCR enrichment
- Single primer PCR
- Inverse PCR
- **Capture + probe extension**
- Paired-end sekvenace
- Germinální - >50x coverage




DNA re-sequencing

Somatic variants

- Tumor x Normal
 - WGS, WES, panels
 - Tumor only variants
- Tumor only
 - Panels, genes
 - Variants germline i somatic
 - Interpretation according to VAF, annotation
- UMI, Duplex-seq adaptors
 - Identification of PCR duplicates
 - Error correction

Levels of error correction and sensitivity





A female scientist with long red hair is wearing a black head-mounted display (HMD) and a blue lab coat. She is examining a green plant with yellow flowers. The background shows a laboratory setting with other people and equipment. Overlaid on the image are several white icons: a DNA double helix, a brain, a cluster of circles, a gear, and an atom symbol.

Thank you
for your attention



CEITEC



@CEITEC_Brno