C8953 NMR structural analysis - seminar

Few Basic Concepts & Vector model

Kateřina Peterková, Aleš Novotný

423977@mail.muni.cz

March 13, 2019

Assign correct value of chemical shift to labelled NMR active atoms¹:

Assign correct value of chemical shift to labelled NMR active atoms¹:

¹http://www.chem.wisc.edu/areas/reich/chem605/< □ > < ₫ > < ≣ > < ≣ > < < <

Diastereotopicity¹ Determine the equivalency of geminal protons

$$\begin{array}{c} C_6H_6 \\ Si \\ O \\ H_6C_6 \end{array}$$

$$\begin{array}{c} C_6H_6 \\ Si \\ O \\ H_3C \end{array}$$

$$\begin{array}{c} C_6H_6 \\ Si \\ O \\ H_{MM} \end{array}$$

$$\begin{array}{c} CI \\ H_{MM} \\ H \\ Br \end{array}$$

Diastereotopicity¹ Determine the equivalency of geminal protons

$$\begin{array}{c|c}
H_a H_a & C_6 H_6 \\
\hline
O & Si \\
O & C_6
\end{array}$$


$$\begin{array}{c|c}
H_6 C_6
\end{array}$$


¹http://www.chem.wisc.edu/areas/reich/chem605/< □ > < ₫ > < 毫 > < ② < ②


Determination of regioisomers - note number of signals, splitting


and δ


300 MHz ¹H NMR spectra in CDCl₃













This work by Dr. James S. Nowick, Professor of Chemistry, University of California, Irvine, is licensed under a Creative Commons Attribution 4.0 International License. Spectra are from Signa-Aldrich (www.sigmaaldrich.com) under fair use.

Determination of regioisomers - note number of signals, splitting


and δ

300 MHz ¹H NMR spectra in CDCl₃






This work by Dr. James S. Nowick, Professor of Chemistry, University of California, Irvine, is licensed under a Creative Commons Attribution 4.0 International License. Spectra are from Sigma-Aldrich (www.sigmaaldrich.com) under fair use.




Draw the estimate of ¹³C NMR spectrum (with and without ¹H decoupling)

Draw the estimate of ¹³C NMR spectrum (with and without ¹H decoupling)




Analysis of simple pulse sequences using vector model

- simple model based on rotation of the vector of bulk magnetization in the plane perpendicular to the vector of magnetic field, direction is determined by the "right-hand rule"
- NMR signal is detectable only as coherent magnetization oscillating in xy plane
- the free precession ω (due to the B₀) of magnetization vector is eliminated by introducing rotating frame $\omega_0 \Rightarrow$ magnetic field of excitation pulses (B₁) is motionless and the individual resonance frequencies differs in so called offset $\Omega_i = \omega_i \omega_0$
- applicability of vector model is rather limited to simple single-quantum experiments without transfer of polarisation




T_1 relaxation

Apply following sequence (inversion recovery) to isolated spin characterized by **a)** $\tau = 2 * T_1$ and **b)** $\tau = 0.2 * T_1$. Draw semi-quantitatively resulting spectrum.



T_1 relaxation

Apply following sequence (inversion recovery) to isolated spin characterized by **a)** $\tau = 2 * T_1$ and **b)** $\tau = 0.2 * T_1$. Draw semi-quantitatively resulting spectrum.




1-1 sequence

Draw the evolution of macroscopic magnetization through the sequence:

90(y) -
$$\tau$$
 - 90(y) - aq

Consider the evolution of an isolated spin due to the chemical shift.

- 1. How does the result differ for the following offsets: $\Omega \tau = 0, \pi/2, \pi$.
- 2. Draw lineshapes of resulting signal assuming the a) y+ b) x+ corresponds to zero phase of receiver (prior phase correction).




1-1 sequence

Draw the evolution of macroscopic magnetization through the sequence:

90(y) -
$$\tau$$
 - 90(y) - aq

Consider the evolution of an isolated spin due to the chemical shift.


- 1. How does the result differ for the following offsets: $\Omega \tau = 0, \pi/2, \pi$.
- 2. Draw lineshapes of resulting signal assuming the a) y+ b) x+ corresponds to zero phase of receiver (prior phase correction).



Heteronuclear spin echo

By using vector diagrams determine the result of attached pulse sequence.


- 1. **Ignore 180 pulse** in hydrogen channel for isolated spin systems **a)** ¹³C-¹H and **b)** ¹³C-¹H₂. Explain the role of CPD block.
- 2. Lets consider **the complete sequence** and isolated spin systems **a)** ¹³C-¹H and **b)** ¹³C-¹H₂.



Heteronuclear spin echo

By using vector diagrams determine the result of attached pulse sequence.


1. **Ignore 180 pulse** in hydrogen channel for isolated spin systems **a)** ¹³C-¹H and **b)** ¹³C-¹H₂. Explain the role of CPD block.



Heteronuclear spin echo

By using vector diagrams determine the result of attached pulse sequence.

2. Lets consider the complete sequence and isolated spin systems a) $^{13}C^{-1}H$ and b) $^{13}C^{-1}H_2$.



Next topic

edited 1D 13C spectra, 2D NMR - homonuclear experiments