


C8953
NMR strukturní analýza
seminář

Heteronuclear correlations


Jan Novotný
176003@is.muni.cz

April 6, 2022


NOESY - Eserine in DMSO




NOESY - Eserine in DMSO




NOESY - Eserine




NOESY - Eserine



1D ^{13}C NMR



100 MHz, CDCl_3



Polarization transfer

- ▶ bigger population difference of ^1H nucleus is transferred via J-coupling to less sensitive nucleus X (^{13}C , ^{15}N)
- ▶ fundamental building block of heteronuclear correlation experiments: in 2D-HX experiment each crosspeak manifests interaction of H and X nucleus coupled through bonds
- ▶ **Task: Draw the evolution of magnetization during basic INEPT pulse sequence. Consider C-H interacting pair.**

HMQC (Heteronuclear Multiple Quantum Correlation)

HSQC (Heteronuclear Single Quantum Correlation)



correlate ^1H -X ($\text{X} = ^{13}\text{C}, ^{15}\text{N}, \dots$) based on
 $^1J_{\text{HX}}$

HMQC (a)

- + more robust experiment
- + change of parameters - HMBC
- lower sensitivity and worse resolution

HSQC (b)

- + better resolution, sensitivity
- + part of more complex multidimensional experiments
- less robust



Practical notes ^1H -X HSQC

- ▶ resolution of overlaps
- ▶ routine experiments to control biomolecular sample
- ▶ easy identification of geminal protons
- ▶ indirect determination of protons bonded to NMR inactive heteroatom
- ▶ heteronuclear correlation \Rightarrow no diagonal crosspeak, no symmetry
- ▶ X decoupled during acquisition \Rightarrow singlet crosspeak

HMBC(Heteronuclear Multiple-Bond Correlation)


heteronuclear correlation based on long-range H-X spin-spin interaction($^nJ_{HX}$, $n>1$)

- ▶ utilizes polarization transfer from H through 2-5 bonds on heteroatom (^{13}C , ^{15}N)
- ▶ allows to detect quaternary heteroatoms (Cq) or connect signals among isolated spin systems


HMBC

correct settings of d1, d2 fo evolution of *J*-coupling necessary


- ▶ $d1=1/2*^1J_{C-H}$ - (120-180 Hz)
- ▶ $d2=1/2*^{2-5}J_{C-H}$ - (3-12 Hz)




^1H - ^{13}C HMBC + ^1H - ^{13}C HSQC




^1H - ^{13}C HMBC + ^1H - ^{13}C HSQC



^1H - ^{13}C HMBC + ^1H - ^{13}C HSQC





^1H - ^{13}C HMBC + ^1H - ^{13}C HSQC





Folic Acid: ^1H 1D

1H FOLATE




Folic Acid: ^{13}C 1D


^{13}C FOLATE




Folic Acid: ^1H - ^{15}N HMBC




Folic Acid: ^1H - ^{15}N HMBC



Folic Acid: ^1H - ^{13}C HMBC



Folic Acid: ^1H - ^{13}C HMBC



Next session:

Complex tasks