

PREPARATION, PHYSICO-CHEMICAL PROPERTIES AND APPLICATION OF QUANTUM DOTS IN BIOANALYSES

Ivona Svobodová

Institute of Analytical Chemistry of the ASCR, v.v.i., Veveří 97, 602 00 Brno, Czech Republic
svobodova@iach.cz


Semiconductor nanocrystals – quantum dots

- semiconductor nano-scale inorganic crystals (1-10 nm)
- core from elements of II. and VI. or III. and V. group
- surface modified by polar organic molecules to increase hydrophilicity



Quantum dots - modification

- different materials (ZnS; ZnSe; PbS; CdS; CdSe; CdTe)
- surface coated (CdSe/ZnS, polymer)
- silanization
- ligands (TOP/TOPO, oleic acid, dithiotreitol, thioglycol, 3-mercaptopropionic acid, 2-mercaptoethylamine)



J. AM. CHEM. SOC. 2006, 128, 10171-10180 9 10171
Hao Zhang, Dayang Wang, Bai Yang

Quantum dots - preparation

1st step: preparation of hydrogen telluride

$$4 \text{ NaBH}_4 + 2 \text{ Te} + 7 \text{ H}_2\text{O} \rightarrow 2 \text{ NaHTe} + \text{Na}_2\text{B}_4\text{O}_7 + 14 \text{ H}_2$$

a) 6 hours 0 °C
b) 30 minute 80 °C
c) over night 2 - 8 °C




2nd step: quantum dots formation

$$\text{CdCl}_2 + \text{NaHTe} + \text{MPA or MA} + \text{heating}$$

MPA 3-mercaptopropionic acid: HS-CH₂-CH₂-COOH
MA 2-mercaptoethylamin: HS-CH₂-CH₂-NH₂

Quantum dots - coating



CdTe/CdS

$$\text{CdTe} + \text{CdCl}_2 + \text{Na}_2\text{S} + \text{MPA or MA} + \text{heating}$$

CdTe/CdS/ZnS

$$\text{CdTe/CdS} + \text{ZnCl}_2 + \text{Na}_2\text{S} + \text{MPA or MA} + \text{heating}$$







Quantum dots - preparation







Dependence of QDs size on reflux time







- different reflux times in second reaction step
- increase of emission wavelength and particle diameter with reflux time















Refluxing time	10 min	1 h	2.5 h	13.5 h	16.5 h	27 h	41 h	44 h
Emission maximum (nm)	498	517	536	652	667	701	737	750
Particle size (nm)	1.99	2.67	2.86	4.02	4.17	4.80	4.87	5











Quantum dots – application

Vizualization of organelles in cells

- analyte migration monitoring in cells, pathogenes detection

HeLa (top left), HeLa (top right), nucleus – 655 nm (pink overlay), Golgi apparatus – 585 nm (yellow overlay), microtubules – 525 nm (green overlay).
 Hep-2 cells (bottom left), Mouse small intestine (bottom right).

Quantum dots – application

Cancer tissues labeling

- imaging and localization of cancer tissues
 - using conjugated QDs with specific cancer antibodies

Digital image (top left), Untagged (top right), Diffuse fluorescence (bottom left), Quantum Dots (bottom right).

Quantum dots – application

Solar cells

- Solar cells from CdSe nanorods (Prof. A.Paul Alivisatos, Lawrence Berkeley National Laboratory)
 - QDs absorbed broader spectrum of sun light
 - QDs allows deposit photosensitive layer on wide range of materials
 - price

Acknowledgment

Foret František, Klepářík Karel, Hezinová Věra, Lišková
 Marcela, Přikryl Jan – Institute of Analytical Chemistry AS CR
HRTEM
 Dr. Mariana Klementová –Institute of Inorganic Chemistry AS CR
Fluorescence and fluorescence lifetimes measurement
 Mgr. Petru Táborskému, Ph.D., doc. Janu Preislerovi, Ph.D. and
 doc. Přemysl Lubal, Ph.D.
 – Department of Chemistry, Faculty of Science, Masaryk University
Financial support
 GAAV - KAN400310651 and KJB400310709
 GACR - GA203/08/1680
 MSMT - LC06023
 AV0Z40310501

