

Document Name: Lenka\_Gabka 24.11.dopo.sds  
Plate Type: Absolute Quantification  
User: 7500RTPCR

Document Information

Operator: zdenek  
Run Date: Wednesday November 24 2010 15:53:31  
Last Modified: Wednesday November 24 2010 17:13:47  
Instrument Type: Applied Biosystems 7500 Real-Time PCR System

Comments:  
SDS v1.3.1

Thermal Cycler Profile

Stage	Repetitions	Temperature
1	1	50.0 °C
2	1	95.0 °C
3	40	95.0 °C 60.0 °C

Standard 7500 Mode

Data Collection : Stage 3 Step 1  
PCR Volume: 25 µL

Well	Sample Name	Detector
A1	NTC L	GAPDH cycD1
A2	NTC	GAPDH cycD1
A3	NTC	GAPDH cycD1
A4	NO RT	GAPDH cycD1
A5	NO RT	GAPDH cycD1
A6	NO RT	GAPDH cycD1
A7	cDNA	GAPDH cycD1
A8	cDNA	GAPDH cycD1
A9	cDNA	GAPDH cycD1
A10	EtOH 0,5	GAPDH cycD1
A11	EtOH 0,5	GAPDH cycD1
A12	pouze MM	GAPDH cycD1
B1	pouze MM	GAPDH

B2	pouze MM	cycD1 GAPDH
B3	pouze MM	cycD1 GAPDH
C1	NTC G	cycD1 GAPDH
C2	NTC	cycD1 GAPDH
C3	NTC	cycD1 GAPDH
C4	NO RT	cycD1 GAPDH
C5	NO RT	cycD1 GAPDH
C6	NO RT	cycD1 GAPDH
C7	cDNA	cycD1 GAPDH
C8	cDNA	cycD1 GAPDH
C9	cDNA	cycD1 GAPDH
C10	pouze MM	cycD1 GAPDH
C11	pouze MM	cycD1 GAPDH
C12	pouze MM	cycD1 GAPDH
D1	EtOH 1	cycD1 GAPDH
D2	EtOH 1	cycD1 GAPDH
D3	pouze MM	cycD1 GAPDH

Time	Ramp Rate	Auto Increment
2:00		100
10:00		100
0:15		100
1:00		100

Task	Ct	StdDev Ct	Qty	Mean Qty	StdDev Qty	Filtered	Tm
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	19.6	0.875					
Unknown	31	0.441					
Unknown	19.61	0.875					
Unknown	30.99	0.441					
Unknown	20.07	0.875					
Unknown	31.35	0.441					
Unknown	31.51	5.307					
Unknown	Undetermined						
Unknown	24	5.307					
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						
Unknown	Undetermined						

