

# CG920 Genomics

## Lesson 4 Forward Genetics

Jan Hejátko

**Functional Genomics and Proteomics of Plants,**  
Mendel Centre for Plant Genomics and Proteomics,  
Central European Institute of Technology (CEITEC), Masaryk University, Brno  
[hejatko@sci.muni.cz](mailto:hejatko@sci.muni.cz), [www.ceitec.muni.cz](http://www.ceitec.muni.cz)



INVESTICE DO ROZVOJE VZDĚLÁVÁNÍ

Tato prezentace je spolufinancována  
Evropským sociálním fondem  
a státním rozpočtem České republiky

# Outline

- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype
    - metabolic profile
    - expression of genes of interest
  - Identification of the mutated locus
    - plasmid rescue
    - iPCR
- Use of libraries of point mutants in forward genetics
  - Positional cloning

# Outline

- Forward vs. reverse genetics



# „Classical“ genetics *versus* „reverse genetics“ approaches in functional genomics

## RANDOM MUTAGENESIS

### „Classical genetics“ approach

EMS



1. IDENTIFICATION OF PHENOTYPE

2. GENE MAPPING

3. GENE IDENTIFICATION  
- position cloning

$h \times n$



(retro)transposons

### „Reverse genetics“ approach

T-DNA



1. ISOLATION OF SEQUENCE-  
-SPECIFIC MUTANT

2. IDENTIFICATION OF  
PHENOTYPE



3. PROOF OF CAUSAL RELATIONSHIP  
BETWEEN INSERTION AND  
PHENOTYPE

# Outline

- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype



# Insertional mutagenesis in forward genetics approaches

- Use of insertional mutagenesis for study of carcinogenesis
  - Infection of E $\mu$ Myc mice by MoMuLV retrovirus leads to lymphomas formation, which arose due to activation of Pim kinases (40 % activation of Pim1, 15 % activation of Pim2), molecular targets of these kinases were unknown





# Insertional mutagenesis in forward genetics approaches

- Use of insertional mutagenesis for study of carcinogenesis
  - Infection of E $\mu$ Myc *pim1* mutants by MoMuLV retrovirus leads to lymphomas formation, which in 90 % contain insertion nearby (activation) *Pim2*



# Insertional mutagenesis in forward genetics approaches


- Use of insertional mutagenesis for study of carcinogenesis
  - Infection of  $E\mu$ Myc double mutants *pim1*, *pim2* by MoMuLV retrovirus leads to lymphomas formation, which can be expected to activate either one of the signalling partner of Pim proteins (Y), one of the proteins of Pim signalling pathway (X) or to activate some of the related pathways leading to lymphomagenesis (Z).






# Insertional mutagenesis in forward genetics approaches

- Isolation of genomic regions adjacent to the insertion site of the provirus
  - Cleavage of genomic DNA and ligation of special linkers, so-called *splinkerettes* (increasing the specificity of amplification)




Devon et al., Nucl Acid Res (1994)




Mikkers et al., Nature Gen (2002)

# Insertional mutagenesis in forward genetics approaches

- Isolation of genomic regions adjacent to the insertion site of the provirus
  - First amplification using specific primers




Devon et al., Nucl Acid Res (1994)




Mikkers et al., Nature Gen (2002)

# Insertional mutagenesis in forward genetics approaches

- Isolation of genomic regions adjacent to the insertion site of the provirus
  - Second amplification using nested primers (increasing the specificity)




Devon et al., Nucl Acid Res (1994)




Mikkers et al., Nature Gen (2002)

# Insertional mutagenesis in forward genetics approaches

- Isolation of genomic regions adjacent to the insertion site of the provirus
  - Sequencing and localization of regions adjacent to provirus by searching in annotated databases of mouse genome



Devon et al., Nucl Acid Res (1994)



Mikkers et al., Nature Gen (2002)

In case of splinkerette, the primer is of the same sequence as the top strand and therefore it is unable to act as a primer until the complement of this strand has been synthesized (from the insert-specific primer at the right-hand side).

# Outline


- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype
    - **metabolic profile**

# Metabolic profiling

- Metabolic profiling of plants
  - Mass and automated analysis of metabolites (up to 25,000) by GC-MS techniques in libraries of T-DNA mutants




Metabolite	RT	MS	MS	MS	MS	MS	MS	MS	MS			
1. alanine	100.1	0.32	175.7	0.77	845.7	833	230.7	0.28	100.7	0.15	531.7	8.36
2. aspartate	100.1	0.39	146.1	0.96	291.7	0.24	188.7	0.28	100.7	0.20	348.7	0.24
3. glutamate	100.1	0.39	131.1	0.94	100.7	0.20	137.7	0.20	100.7	0.20	12.7	0.07
4. aspartate	100.1	0.37	128.7	0.86	100.7	0.20	104.7	0.20	100.7	0.20	1.7	0.01
5. valerate	100.1	0.38	304.7	0.20	137.7	0.13	131.7	0.10	100.7	0.20	1.7	0.01
6. succinate	100.1	0.38	83.7	0.88	335.7	0.81	137.7	0.22	100.7	0.20	1.7	0.01
7. lactate	100.1	0.36	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.01
8. pyruvate	100.1	0.37	122.7	0.37	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.01
9. glutamate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.01
10. glycine	100.1	0.31	95.7	0.37	95.7	0.13	95.7	0.13	100.7	0.20	1.7	0.01
11. malate	100.1	0.31	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.01
12. homocysteine	100.1	0.31	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.01
13. histidine	100.1	0.32	84.7	0.88	397.7	0.28	131.7	0.18	100.7	0.20	6.6	0.48
14. valerate	100.1	0.39	100.7	0.20	93.7	0.14	127.7	0.28	100.7	0.20	6.2	0.12
15. malate	100.1	0.32	100.7	0.20	111.7	0.19	103.7	0.08	100.7	0.20	1.7	0.02
17. citrate	100.1	0.38	304.7	0.81	845.7	0.88	330.7	0.28	100.7	0.20	6.7	0.14
18. malonate	100.1	0.33	104.7	0.33	344.7	0.11	236.7	0.11	100.7	0.20	6.6	0.13
19. succinate	100.1	0.32	100.7	0.20	845.7	0.88	841.7	0.28	100.7	0.20	6.7	0.14
20. succinate	100.1	0.34	100.7	0.20	837.7	0.86	840.7	0.08	100.7	0.20	6.6	0.18
21. succinate	100.1	0.38	100.7	0.20	111.7	0.21	100.7	0.20	100.7	0.20	1.7	0.02
22. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
23. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
24. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
25. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
26. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
27. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
28. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
29. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
30. succinate	100.1	0.38	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
1. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
2. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
3. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
4. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
5. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
6. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
7. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
8. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
9. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
10. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
11. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
12. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
13. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
14. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
15. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
16. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
17. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
18. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
19. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
20. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
21. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
22. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
23. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
24. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
25. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
26. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
27. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
28. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
29. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02
30. succinate	100.1	0.32	100.7	0.20	100.7	0.20	100.7	0.20	100.7	0.20	1.7	0.02



Continued

# Metabolic profiling

- Metabolic profiling of plants
  - Mass and automated analysis of metabolites (up to 25.000) by GC-MS techniques in libraries of T-DNA mutants
  - Identification of interesting (even commercially interesting) mutants



# Metabolic profiling


- Metabolic profiling of plants
  - Mass and automated analysis of metabolites (up to 25.000) by GC-MS techniques in libraries of T-DNA mutants
  - Identification of interesting (even commercially interesting) mutants
  - Fast and easy isolation of genes through identification of sequences affected by T-DNA





# Metabolic profiling

- Metabolic profiling of plants
  - Possibility to use special techniques, e.g. microdissection





# Outline

- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype
    - metabolic profile
    - **expression of genes of interest**

# Expression profile

- Identification of mutants with a change in the expression profile
  - Analysis of expression profile (pattern) of the gene and identification of mutants with altered pattern of expression




# Expression profile

- Identification of mutants with a change in the expression profile
  - Analysis of expression profile (pattern) of the gene and identification of mutants with altered pattern of expression
  - Possibility of partial automation (virtual digital microscopy)



# Expression profile




# Outline

- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype
    - metabolic profile
    - expression of genes of interest
  - Identification of the mutated locus
    - plasmid rescue
    - iPCR

# Identification of mutated locus

- Identification of chromosomal rearrangements responsible for bushy phenotype of *Arabidopsis*
  - Description of phenotype


# Identification of mutant



- Crinkled leaves
- Bushy phenotype (branching defective)
- No trichomes on leaves and stems
- Late senescence



# Identification of mutant




- Male sterility, defects in stamen filament elongation (A,B) (compare with wild type C)

# Identification of mutated locus


- Identification of chromosomal rearrangements responsible for bushy phenotype of *Arabidopsis*
  - Description of phenotype
  - Identification of T-DNA mutated region

# Identification of mutated locus

## 1. Identification of region of genomic DNA adjacent to the *left border* using *plasmid rescue*





- Restriction digestion (*EcoRI*) of mutant genomic DNA
- Religation and transformation of *E. coli*
- Isolation of plasmid DNA from positively selected clones
- Identified sequence was identical to gene for *NAD7* coded by mtDNA



# Identification of mutated locus

## 2. Identification of region of genomic DNA adjacent to the **right border** using **inversion PCR (iPCR)**




# Identification of mutated locus

- Identification of chromosomal rearrangements responsible for bushy phenotype of *Arabidopsis*
  - Description of phenotype
  - Identification of T-DNA mutated region
  - Localization of T-DNA insertion site in *Arabidopsis* genome

# Searching in library IGF-BAC

- Genome library containing 10.752 clones with an average size of an insert of 100 kb
- Bacterial clones arranged in the microtiter plates
- Library loaded onto nylon filters for hybridization with the radiolabeled probe



# Mapping with IGF-BAC database

## I. Sequences adjacent to the left border of T-DNA


- 28 positively hybridizing clones in total
- 19 of them located on chromosome 2
- 18 of them similar with mtDNA

## II. Sequences adjacent to the right border of T-DNA


- 6 positively hybridizing clones in total
- all of them located on chromosome 2

# Localization of genomic T-DNA adjacent to both left and right T-DNA borders on chromosome 2

## Sequences adjacent to *right* and *left* border of T-DNA



- There was probably an inversion of almost entire chromosome 2





# Outline

- Forward vs. reverse genetics
- Use of libraries of insertional mutants in forward genetics
  - Searching in libraries of insertional mutants according to:
    - anatomically or morphologically detectable phenotype
    - metabolic profile
    - expression of genes of interest
  - Identification of the mutated locus
    - plasmid rescue
    - iPCR
- Use of libraries of point mutants in forward genetics
  - Positional cloning







# Identification of mutated locus

## ■ Positional cloning

- Principle: co-segregation analysis of segregating population (mostly of offspring of backcrosses) with molecular markers
- **SSLP** (Simple Sequence Length Polymorphism)
  - Polymorphism of genome (PCR products) length, amplified using specific primers
- **RFLP** (Restriction Fragment Length Polymorphism)
  - Detection by Southern blot (PCR after digestion of the genomic DNA and ligation of adapters)
- **CAPS** (Cleaved Amplified Polymorphic Sequence)
  - Restriction fragment length polymorphism, genome segments amplified by PCR
- **RAPD** (Randomly Amplified Polymorphic DNA)
  - Polymorphism of length of randomly amplified genome segments, using short 8-10bp primers

# Positional cloning

Preparation of mapping population




**Recombinant analysis** – determining the percentage of recombination between mutation and molecular marker

$$r [\%] = \frac{\text{number of chromosomes of Col}}{\text{number of all the chromosomes}} \times 100$$

**F2 mutants**

Ler Col




marker I – linked

5 mutants

$$1/10 \times 100 = 10\%$$

**F2 mutants**

Ler Col




marker II - no linkage

6 mutants

$$7/12 \times 100 = 58\%$$

- Analysis of approximately 2000 mutant lines
- Determining the closest (still) segregating marker
- Identification of mutation by sequencing

# Map of DNA molecular markers



# Markers for fine mapping

- AGI Map
- Lister & Dean RI
- Classical
- mi-RFLP
- Goodman
- GoodmanBAC
- TIGR
- Finkelstein
- Altmann

## Maps for Chromosome 2

for all Maps: [Search Options:](#)

Selected Maps   
     
 Display All Rows



[MapViewer Home](#)

[Release Note](#)

[View Print-Version](#)

### AGI Map



[Zoom to:](#)


Zoom up to 200x to see genes!

Search by name (e.g. UFO)

Select range (e.g. 1500-2000)

[AGI Map color key](#)

### Lister & Dean RI



[Zoom to:](#)

Search by name (e.g. UFO)

