

# How to visualize genes and their products

CG920 Genomics  
Lecture 9


Kamil Růžička  
FGP CEITEC MU

# Outline

- reporter genes
- promoter fusions
- visualizing proteins
- visualizing RNA
- dynamics of protein imaging: FRAP, photoactivable proteins, FLIM, FCS


# WHERE YOU SIT IN CLASS/SEMINAR

And what it says about you:



WWW.PHDCOMICS.COM

Proximity to Lecturer:



$$X = \frac{\text{How much you care}}{\text{How sleepy you are}}$$

JORGE CHAM © 2008

# Promoter activity monitoring



1-10 kb prior to ATG

LacZ, GUS

Luciferase

GFP

# Reporter genes

- LacZ, GUS
- Luciferase
- GFP


some need external substrate, some do not

# LacZ, GUS – rhapsody in blue

promoter

LacZ


terminator



(in case of GUS – X-Gluc)

# LacZ, GUS

LacZ/ GUS:




worm, mouse – LacZ, plants – GUS

# Luciferase

promoter

luciferase

terminator





(similar to chemiluminescence)


What's the difference between fluorescence and luminiscence? 8




# Luciferase




# How does fluorescence work?




# How does a fluorescence microscope work?




# Stokes shift



George G. Stokes



# How does a confocal microscope work?



What are advantages of confocal microscopy?

# Live imaging

GFP discovery - Nobel Prize 2008




Osamu Shimomura



Martin Chalfie

Roger Tsien

# Fluorescent proteins on the market (Tsien's fruits)




# Excitation and emission





# Multicolored fluorescent protein image (neurons)




# Promoter-GFP

promoter

GFP

terminator



# Promoter activity monitoring

choice of suitable reporter

- LacZ, GUS
- Luciferase
- GFP

accessibility, sensitivity, accuracy...

# Promoter activity monitoring


- LacZ, GUS
  - easy assay, also on sections, easy imaging
  - substrate must diffuse, kills the organism
- luciferase
  - good quantification, very sensitive, no autofluorescence
  - substrate must diffuse, special machine, dark
- GFP
  - good sensitivity, colocalization with other dyes/promoters possible, no substrate needed
  - only in vivo, autofluorescence, thin transparent sample; it should be ER localized in plants

# Luminiscent mouse better than fluorescent mouse

## *In Vivo* Comparison of Bioluminescence and Fluorescence (I.M.)

- Fluorescent signal is limited by tissue autofluorescence
- The bioluminescent signal level is ~300x lower, yet the signal to background is 160x higher

### Bioluminescence




Background flux  $\sim 2.6 \times 10^3$  p/s  
Signal flux  $\sim 2.8 \times 10^6$  p/s  
Signal/background  $\sim 1100$   
Min. detectable cells  $\sim 900$

Background flux  $\sim 1.2 \times 10^8$  p/s  
Signal flux  $\sim 8.3 \times 10^8$  p/s  
Signal/background  $\sim 6.7$   
Min. detectable cells 150,000

Left:  $1 \times 10^6$  HeLa-luc/PKH26 cells  
Right:  $1 \times 10^6$  HeLa-luc cells

### Fluorescence



# Promoter activity monitoring

## Pros:


- easy to clone, easy to visualize
- usually some signal seen – cheers you up!
- can be used in less accessible organs

## Cons:


- limited information about gene product (mRNA, protein etc.)
- needs cloning and transformation
- neglects regulatory elements (introns, UTRs etc.)
- length of promoter given arbitrarily

# Translational GFP fusions

N-terminal fusion




C-terminal fusion



fusion inside the coding sequence





# GFP and membrane proteins




It is good to have GFP tag localized inside the cell (plants)




# Expression of isoforms



# Isn't this better?



# Expression of isoforms



# Fluorescent protein fusion

## Pros:

- in vivo imaging


## Cons:

- not always functional
- transformation needed
- transparent material (you can sometimes fix GFP signal, however)
- sometimes GFP artifacts (tag doesn't allow proper targeting)

# Why to visualize all this stuff



promoter



translational


# Why to visualize all this stuff



*pSHR :: GFP*

promoter


- 1 – epidermis
- 2 – cortex
- 3 – endodermis
- 4 – stele



*pSHR::SHR::GFP*


translational

# Why to visualize all this stuff



**BANG! SHR moves from stele to endodermis**

# Protein immunolocalization



Most favorite animals:



- rabbit (too many rabbits)
- mouse (low volume)
- goat
- chicken
- rat
- sheep
- donkey
- guinea pig

secondary: anti-rabbit from no-rabbit, anti-mouse from no-mouse, etc.




# Protein immunolocalization

## immunolocalization - fluorescently



# Protein immunolocalization



immunolocalization



Fluorescent dyes conjugated to secondary (examples):

- FITC (obsolete)
- CY3, CY5
- Alexa (488, 568, 633)

# Fluorescent Dyes and Proteins



Dye Name	Excitation Max	Emission Max

# Protein immunolocalization


## Pros:

- no need to clone or transform or cross
- direct (if no tag used)
- allows sectioning (less accessible tissues)


## Cons:

- fixed material only
- excellent antibodies only, sometimes tricky

# GFP tag partially retains PIN1 in endoplasmic reticulum (-> artifact)





PIN1-GFP



anti-PIN1



# Protein localization - immunogold

immunolocalization - immunogold



electron microscope

# Immunogold colocalization



# Pros/cons

## Pros:

- direct
- nothing can beat the resolution



## Cons:

- very tricky (needs rather expert)
- huge experience for interpretation needed




Can we visualize posttranslational modifications?

# Also RNA can be visualized



# Localization of mRNA

## RNA hybridization *in situ*



# Visualization of mRNA

## RNA hybridization *in situ*

### Pros

- classical technique in developmental biology
- no transgenes needed

### Cons

- tedious, tricky, no success guaranteed
- only on fixed samples



For shorter RNAs (miRNA etc.):

- LNA probes needed




# Single-molecule detection using multiple probes

- ~ 48 oligonucleotide probes provide sufficient signal to detect a single mRNA molecule



# Also mRNA can be visualized *in vivo*




Ash1 mRNA localized to the tip of the daughter cell

# $\lambda N_{22}$ system – RNA imaging *in vivo*

nuclear localization signal


viral RNA binding protein



# Drawbacks of $\lambda N_{22}$ system - we have SPINACH

GACGCAACUGAAUGAAA  
 UGGUGAAGGACGGGUCC  
 AGGUGUGGCUGCUUCGG  
 CAGUGCAGCUUGUUGAG  
 UAGAGUGUGAGCUCCGU  
 AACUAGUCGCGUC


+



RNA fusion



aptamer





blue-DNA


green-RNA



# Other vegetables than SPINACH



# Advanced confocal techniques




# (slightly) Advanced confocal techniques

- FRAP
- photoactivatable FP
- FCS



# FRAP

## Fluorescence Recovery After Photobleaching


region of interest (ROI)



# FRAP




you can quantify fluorescence..  
(ImageJ is our friend)






	mean	min	max
A	90.404	49	113
C	8.556	3	8
D	39.934	19	63





# FRAP – bleaching curve



What does the curve tell?





# FRAP – bleaching curve




# iFRAP

inverse FRAP






# iFRAP – dissociation of premRNA from speckles



# FRAP derivatives

## FLIP

Fluorescence Loss After Photobleaching




- bleaching process is repeated during the experiment
- for studying general protein turnovers in compartments
- scientific question here: is there a fraction of protein which does not leave the bright green patches

# FRAP derivatives

## FLAP

Fluorescence Localization after Photobleaching




- two fluorochromes on one protein– one bleached, non bleached as control

# Perhaps better scheme than previous

CFP not bleached

YFP bleached



RED=CFP-YFP


# FRAP - advantages

- not only proteins (also other dyes)
- tells you more than simple life imaging movie

# FRAP – pitfalls

- your cells are moving
- high energy needed to bleach the ROI
  - long time needed to bleach
  - can damage your material
- usually only one ROI can be observed – time consuming
- for gourmets perhaps awkward (although more reliable and robust)








# Photoactivatable fluorescent proteins



photoactivation  
(UV)

aquaporin PIP2  
undergoes  
lateral diffusion

# Photoactivatable proteins

	Protein		Cell or organelle				
	Repeated short-term tracking	Long-term tracking	In culture		In tissue		
	<b>Dronpa</b>	<b>PAmRFP1</b>	<b>mEosFP</b>	<b>PA-GFP</b>	<b>PS-CFP2</b>	<b>Kaede</b>	<b>KFP1</b>
	Monomer	Monomer	Monomer	Monomer	Monomer	Tetramer	Tetramer
	Reversible	Irreversible	Irreversible	Irreversible	Irreversible	Irreversible	(Ir)reversible
Fluorescence changes during photoactivation							
High brightness	✓			✓		✓	
High contrast			✓		✓	✓	
Dual labelling with red and green fluorescent proteins	✓	✓		✓	✓		✓
Low phototoxicity of the activation light							✓

Copyright © (2005) Nature Publishing Group  
 Nature Reviews | **Molecular Cell Biology**

Dronpa, Kaede, Eos – probably the most popular



# Photoactivable proteins

## Advantages:

- elegant, can be convincing

## Disadvantages:

- very weak signal
- each material needs optimization

# Remarks

- your material is 3D
- protein *de novo* synthesis in some experiments (e.g. cycloheximide stops translation)


# FLIM

Fluorescence Life Time Imaging Microscopy

## Fluorochromes

- excitation spectra
- emission spectra
- **unique lifetime**

# FLIM - applications




# FLIM - applications

Lifetime sensitive to almost everything:

- pH
- ionic strength
- solution polarity
- other fluorochrome


Protein-protein interactions  
(FRET-FLIM)

# FLIM



indeed, salt changes fluorophore life time  
(American cockroach glands)

# FLIM - discrimination of autofluorescence



(be careful with the interpretation)

Q: What might be the easier experiment to confirm autofluorescence?


# FLIM

- experience needed
- special module on your confocal needed



# FCS



## Fluorescence Correlation Spectroscopy



It is counted, how many times the fluorescent molecule comes through the focal plane.

Autocorrelation analysis: the way how to discriminate the diffusions speeds of particles.



# FCS



# FCS (FCCS)

fluorescence cross-correlation spectroscopy


$$G_{CC}(\tau) = \frac{\langle I_A(t) \cdot I_B(t + \tau) \rangle}{\langle I_A(t) \rangle \langle I_B(t) \rangle} - 1$$



# Literature

- Paige et al., RNA Mimics of Green Fluorescent Protein, Science 333, 642-646, 2011 (SPINACH and other vegetables)
- <https://www.micro-shop.zeiss.com/index.php?s=452278238ae52c&l=en&p=de&f=f&a=d> (comprehensive and broad list of phluorochromes)
- <http://www.illuminatedcell.com/> - nicely done pages about plant cell imaging
- Ishikawa-Ankerhold et. al. Advanced Fluorescence Microscopy Techniques — FRAP, FLIP, FLAP, FRET and FLIM, Molecules 2012, 17, 4047-4132
- Sambrook & Russell Molecular Cloning: A Laboratory Manual, Third Edition, Cold Spring Harbor Laboratory Press; 3rd edition (January 15, 2001), 13.5 pounds weight

# Photon bunching, if someone asks



Photon detections as a function of time for a) antibunched, b) random, and c) bunched light