

Isothermal titration calorimetry

S2004 – METHODS FOR CHARACTERIZATION OF BIOMOLECULAR
INTERACTIONS: CLASSICAL VERSUS MODERN

Mgr. MONIKA KUBÍČKOVÁ, Ph.D.

SLIDES BY: Mgr. JITKA HOLKOVÁ



MUNI

Outline:

- ▶ Historical background
- ▶ Theory
- ▶ Study of interactions
- ▶ Instrumentation
- ▶ ITC data analysis

S 2004 Methods for characterization of biomolecular interactions

CF BIC – Masaryk University

History of calorimetry

- ▶ Calorimetry – Latin *calor* – heat, Greek *μέτρον* (-metry) – to measure
 - thermodynamic technique based on measurement of heat that may be generated (**exothermic process**) or consumed (**endothermic process**) by sample
- ▶ Calorimeter – instrument for measuring the quantity of heat **released** or **absorbed** in process of chemical reaction

Calorimetric units

- ▶ A single calorie is the amount of energy required to increase the temperature of 1g of water by 1°C.
- ▶ A single joule is the amount of energy required to apply a force of 1 Newton over one meter of distance.
- ▶ 1 calorie = 4.184 J
1 Calorie = 1 kcal = 4184 J
1 J = 0.000239 kcal = 0.2390 cal


History of calorimetry: „Founding Fathers“

▶ **Joseph Black (1728 – 1799)**

- „founder of the calorimetry“
- first scientist who recognized the distinction between heat and temperature

▶ **Antoine Lavoisier (1743-1794)**


▶ **Pierre- Simon Laplace (1749-1827)**



S 2004 Methods for characterization

CF BIC – Masaryk University

First calorimeter - small guiney pig inside





S 2004 Methods for

interactions
Masaryk University

Calorimetry

- ▶ **INDIRECT CALORIMETRY –**
calculates the heat generated by living organism when their metabolic processes yield waste carbon dioxide
- ▶ **DIRECT CALORIMETRY –**
measures heat generated by living organism by placing the entire organism inside the calorimeter for the measurement



Microcalorimetry in cube:


Microcalorimetry

- Direct measurement of heat change (ITC)
- Direct measurement of melting transition temperature to predict thermal stability (DSC)




Broad dynamic range

- Native molecules in solution (biological relevance)
- Very sensitive to accommodate range of affinities



Information rich

- All binding parameters (affinity, stoichiometry, enthalpy and entropy) in a single ITC experiment



Ease-of-use

- No labeling or immobilization necessary
- Wide range of solvent/buffer conditions




S 2004 Meth

actions
University

Microcalorimetry


▶ Differential scanning calorimetry DSC

- Biomolecular stability in solution
- Provides insights into mechanisms of unfolding and refolding
- Midpoint (T_m) determination



▶ Isothermal titration calorimetry ITC



- Heat is released or absorbed as a result of the redistribution and formation of non-covalent bonds when the interacting molecules go from the free to the bound state.



With isothermal titration calorimetry you can...

- ▶ Get quick K_D s
- ▶ Measure target activity (Stoichiometry, active concentration)
- ▶ Protein batch activity comparison
- ▶ Confirm drug binding to target
- ▶ Use thermodynamics to guide lead optimization
- ▶ Measure enzyme kinetics

How does it work?



The DP is a measured power differential between the reference and sample cells to maintain a zero temperature between the cells




DP = Differential power
 ΔT = Temperature difference

S 2004

CF BIC - Masaryk University

Basics of ITC experiment

Universal technique based on heat detection



Integration of heats are used to extract affinity (K_D), stoichiometry (N) and binding enthalpy (ΔH) using appropriate binding model

The energetics

$$\Delta G = RT \ln K_D$$

$$\Delta G = \Delta H - T\Delta S$$

ΔH , **enthalpy** is indication of changes in hydrogen and van der Waals bonding

$-T\Delta S$, **entropy** is indication of changes in hydrophobic interaction and/or conformational changes

ΔG = Gibbs free energy


ΔH = Enthalpy

ΔS = Entropy

R = Gas constant = $1.985 \text{ cal K}^{-1} \text{ mol}^{-1}$


T = Temperature in Kelvin = $273.15 + t \text{ } ^\circ\text{C}$

K_D = Affinity




The energetics

- ▶ The same affinity and stoichiometry but different enthalpy
- ▶ This tells us there are different binding mechanisms



Same affinity, different energetics!
All three interactions have the same binding energy (ΔG)

- A. Good hydrogen bonding with unfavorable conformational change
- B. Binding dominated by hydrophobic interaction
- C. Favorable hydrogen bonds and hydrophobic interaction




ITC results are used to get insights into mechanism of binding



ITC experiment

Standard set-up:

- ▶ “Ligand” in syringe
- ▶ “Macromolecule” in sample cell
- Reverse arrangement possible
- Concentration and other parameters necessary to set-up the experiment



Performing an ITC experiment




S 2004 Me...

CF D¹-

interactions
aryk University

Competition titration:

Very high and very low affinity systems can be studied using competition titrations




- High affinity ligand added to a solution of the low affinity complex
- High affinity ligand displaces the low affinity ligand
- Change in the apparent affinity and apparent enthalpy
- If parameters for one ligand are known, possible to calculate for the other ligand

S 2004 M€

eractions
aryk University


Single injection method

- ▶ This setup consist of one injection of substrate into enzyme with the aim to reach V_{max} as fast as possible and observe the signal decay associated with substrate depletion and product formation.




Assessment of protein quality by MicroCal™ iTC₂₀₀ system

Peptide binding to protein Batch #1



- ▶ **100%** of Batch #1 protein active based on stoichiometry

Peptide binding to protein Batch #2



- ▶ **23%** of Batch #2 protein active based on stoichiometry

Sample preparation

Sample preparation – „c value“


$$C = [\text{Protein}]/K_D$$

C = 10-100 Great

C = 5-10 and 100-500 Good

C = 1-5 and 500-1000 Accuracy limited

C = < 1 and > 1000 Competitive ITC experiment



S 2004 M

CF BIC


actions
k University

Sample preparation

- ▶ Stability of interacting molecules in conditions of the ITC experiment
- ▶ The cell and syringe buffers must be carefully matched. This is best accomplished by dialyzing both the macromolecule and the ligand in the same buffer.
- ▶ If the ligand is too small for dialysis then dialyze the macromolecule and then dissolve the ligand in the dialyze buffer
- ▶ Accurately measure protein concentration using A280nm

Poor sample preparation = poor data

- ▶ The data show possible difference of measurement of the sample before and after dialysis
- ▶ The large peaks were due to differences in the NaCl concentration between buffers



S 2004 Methods for characterization

CF BIC – Masaryk University

Instrumentation

Malvern Instruments

► Differential scanning calorimetry

MicroCal™ VP-DSC



MicroCal VP-Capillary DSC



► Isothermal titration calorimetry

MicroCal VP-ITC



MicroCal iTC₂₀₀



MicroCal PEAQ™ ITC



MicroCal PEAQ ITC Automated



S 2004 Methods .

- Masaryk University



TA Instruments

▶ Differential scanning calorimetry



▶ Isothermal titration calorimetry

	Standard Volume	Low Volume
Cell Geometry	Fixed Cylindrical	Fixed Cylindrical
Cell Composition	24K Gold / Hastelloy	24K Gold
Active Cell Volume	1.0 mL	190 μ L




ITC – theory of data analysis

What we are going to discuss?

- ▶ Analysis of raw ITC data
- ▶ Data fitting
 - - choice of model
 - - fitting procedure
 - - assessment of goodness of fit
- ▶ Global fit: fitting of multiple data sets



Raw ITC data




✓ **Make a sanity check of the raw data**

Bad quality raw ITC data

But **Buffer mismatch – no dialysis**



Too low reference power



Too short spacing between the injections

S 2004 M

erization

molecular interactions


yk University



Origin-based ITC data analysis software

First steps

- ▶ Adjust integrations
- ▶ Check concentrations
- ▶ Subtract control heats



S 2004 Methods for character

eractions

aryk University


Choice of model

► Models supported by Origin:


► **4 binding models**

- one set of sites
- two sets of sites
- sequential binding
- competitive binding

► **a dimer dissociation model**




“One set of sites” model: parameters defined by binding isotherm




Manual fit initialization : your educated guess

If automatic initialization is not satisfactory, in NL curve fitting box insert your "best-guess" values for parameters and click Chi-Sqr button. A simulated curve will appear next to your experimental data curve. Compare and decide.



The simpler the model, the higher the chance automatic initialization will work.



CF BIC



University

Control of fitting procedure


Tolerance—compare chi-sqr values between two successive iterations

Delta – controls the way partial derivatives are calculated

Parameter constrains – allow to exclude unphysical values of parameters (USE WITH CAUTION)



Quality of fit: dependency of parameters



✓ Dependency value very close to 1 indicates strong cross-correlation and over-parameterization.

S 2004 M

BIC - Masaryk University

Quality of the fit

- Chi-sq
- Parameter dependence
- Errors in the fitted parameters
- Agreement between repeated experiments
- Biochemical and experimental relevance in the parameters returned by the fit

Quality of the fit: fitted parameters

N , number of binding sites

$$Q = \frac{nM_t \Delta H V_o}{2} \left[1 + \frac{X_t}{nM_t} + \frac{1}{nKM_t} - \sqrt{\left(1 + \frac{X_t}{nM_t} + \frac{1}{nKM_t} \right)^2 - \frac{4X_t}{nM_t}} \right]$$


- ▶ “N” is the average number of binding sites per mole of protein in solution, assuming:
 - that all binding sites are identical and independent
 - that you have pure protein (and ligand)
 - that you have given the correct protein and ligand concentrations
 - that all your protein is correctly folded and active

Goodness of the fit: fitted parameters


N, number of binding sites

- ▶ If $N \neq 1$
 - different number of binding sites
 - inaccurate input values for protein and/or ligand concentration
 - protein instability issues
 - compound solubility issues
 - binding does not fit simple independent model


Stoichiometry: Incorrect [Ligand]



N 0.8218
K 6.899E4
 ΔH -1.694E4




N 1.279
K 4.434E4
 ΔH -1.089E4




N 1.023
K 5.543E4
 ΔH -1.361E4


Stoichiometry: Incorrect [Protein]



N 0.8181
K 5.543E4
 ΔH -1.361E4



N 1.278
K 5.543E4
 ΔH -1.361E4






N 1.023
K 5.543E4
 ΔH -1.361E4



- ▶ **Error in syringe concentration results in error in DH, K and N**
- ▶ **Error in cell concentration results in error in N**
- ▶ Put the sample of which you have most control over in the syringe and evaluate accordingly



Reporting results: final figure and parameter box



S 2004

ations
 University


Global fit of multiple datasets collected at different experimental conditions



S 2004 Methou


CT V.

ctions
University



MicroCal PEAQ software-based ITC data analysis

MicroCal PEAQ-ITC software



S 2004 Metho

eractions

ryk University

MicroCal PEAQ-ITC software

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Adjust Fit

Start Analysis Overview Adjust Baseline Assign Controls Adjust Fit Presentation

Experiments Sort by Bin

Experiment Information

Filename	
Temperature (°C)	25.1
[Syr] (M)	7.38e-3
[Cell] (M)	499e-6
Ref. Power (μcal/s)	10.0 (10.4)
Comment	
Control Type	Fitted Offset
Control Experiment	N/A
Results	
Bin Comment	Binding
Competitive Model	No
[Syr] (M)	7.38e-3
[Cell] (M)	499e-6
Ligand In Cell	No
N (sites)	0.771 ± 1.7e-3
K _D (M)	13.9e-6 ± 335e-9
ΔH (kcal/mol)	-5.83 ± 2.1e-2
ΔG (kcal/mol)	-6.63
-TΔS (kcal/mol)	-0.800
Offset (kcal/mol)	-0.157 ± 6.6e-3
Reduced Chi-Sqr (kcal/mol) ²	3.5e-4

Fitting Parameters

Parameter	Vary	Initial Value	Lower Bound	Upper Bound
N (sites)	<input checked="" type="checkbox"/>	0.771	1.0e-3	10.0
K _D (M)	<input checked="" type="checkbox"/>	13.9e-6	1.00e-12	1.00
ΔH (kcal/mol)	<input checked="" type="checkbox"/>	-5.83	-80.0	80.0
Offset (kcal/mol)	<input checked="" type="checkbox"/>	-0.157	-80.0	80.0
[Syr] (M)	<input type="checkbox"/>	7.38e-3	0	1.00
[Cell] (M)	<input type="checkbox"/>	499e-6	0	1.00

Use Competitive Model

Unknown Binder: Strong

[Weak] (M): 499e-6

Known Weak Parameters: Enter Manually

N (sites): 1.00

K_D (M): 1.00e-6


ΔH (kcal/mol): -1.00

Reset Initialize Fit Fit Iterate Once Simplex Fit Show Excluded Injections Show Legend

S 2004

ions
University

MicroCal PEAQ-ITC software



S 2004 Methods for chara

CF BIC – Masaryk University

molecular interactions

MicroCal PEAQ-ITC software – Design experiment

Malvern MicroCal PEAQ-ITC Analysis Software

Analyze Experiment(s) Design Experiment

Mode Guided **Advanced**

Select Experiment Model
One Set of Sites

Experimental Parameters
Temperature* (°C) 25
Reference Power (μcal/s) 10.0
FeedBack High
Stir Speed (rpm)* 750
Initial Delay (s) 60
*not simulated, but saved in method

Injection Details
of Injections 20

Injection	Volume (μL)	Duration (s)	Spacing (s)
1	0.4	0.8	150
2	3.0	6.0	150
3	3.0	6.0	150
4	3.0	6.0	150
5	3.0	6.0	150
6	3.0	6.0	150
7	3.0	6.0	150
8	3.0	6.0	150

Apply to All Apply to Rest

Display Normalized Heat Raw Heat

Concentrations
[Syr] (M) 400e-6
[Cell] (M) 40e-6

Misc.
c-value 40.0 Lock
Total Heat (μcal) 22.3
[Syr]/[Cell] 10.0 Lock

Binding Parameters
N (sites) 1.00
K_D (M) 1.00e-6
ΔH (kcal/mol) -3.00

NOTE: Simulated injection heat variability may not be an accurate representation.

Reset Save As Method

S 2004 M

actions
k University

Thank you for you attention...😊


ITC and **DSC** techniques available in **CF BIC** (Core Facility of Biomolecular Interaction and Crystalization) **bic.ceitec.cz**

CF Head: Prof. Michaela Wimmerová

Contact: Monika Kubíčková, monika.kubickova@ceitec.cz C04/339

Data fit: non-linear least-squares minimisation

- ▶ Fitting procedure evaluates the deviation of the fitted function from the experimental data in terms of chi-squared.
- ▶ In Origin ITC data analysis module minimization is performed iteratively by Marquardt-Levenberg or simplex algorithms



© 2004 Methods for characterization

Interactions
York University

CTV