

GENEROVÁNÍ TĚKAVÝCH SLOUČENIN PRO STOPOVOU PRVKOVOU A SPECIAČNÍ ANALÝZU: VÝHODY A OMEZENÍ

Ústav analytické chemie AVČR, v. v. i.

Oddělení stopové prvkové analýzy



Akademie věd
České republiky


iac
brno

24.11. 2022

Jan Kratzer

GENERATION OF VOLATILE SPECIES

- volatile compounds generation (VCG) / volatile species (VSG)
- **selective conversion of analyte to volatile compound**
- VSG compatible with all spectrometric techniques
- VSG independent of the detector used
- most commonly used for **hydride forming elements (HG)**



As
Sb
Bi
Se
Te
Pb
Sn
Ge

GENERATION OF VOLATILE SPECIES

- **selective conversion of analyte to volatile species**
 - analyte separation from the matrix (interferences minimized)
 - high transport efficiency of analyte into the detector → lower LOD
- analyte **preconcentration** in gaseous phase → further LOD improvement
- **speciation analysis** without chromatography feasible

APPROACHES TO VSG


- **Chemical generation (CVG)**

- analyte reduction by chemical reaction (HCl/NaBH_4)
- high generation efficiency ($\sim 100\%$)



- **Electrochemical generation (EcVG)**

- analyte reduction by current
- low generation efficiency
- potential to reach low LOD



- **Photochemical generation (PVG)**

- analyte reduction by UV radiation



- **Plasma mediated vapor generation (PMVG)**


- interaction with plasma - radicals, excited/metastable species, ions




OUTLINE

- **VSG of non-hydride forming element – Cd**
- **Novel hydride atomizers for AAS based on DBD plasma**
- **VSG-based speciation analysis (Hg, Te, Ge, As)**
- **VSG of mercury based on PMVG of Hg**


- **VSG of non-hydride forming element – Cd**
- Novel hydride atomizers for AAS based on DBD plasma
- VSG-based speciation analysis
- VSG of mercury based on PMVG





conventional approach: HG

2-channel system (HCl/NaBH₄)



VSG of Cd

4-channel system

0.1 mol L⁻¹ KCN
0.001 mol L⁻¹ Cr³⁺

- generation efficiency quantified (**with additives**)
 - by $^{115\text{m}}\text{Cd}$ radioactive indicator $66 \pm 4 \%$
 - by comparison of PN-ICP-MS and VSG-ICP-MS $55 \pm 2 \%$
- HCl/NaBH₄ **without additives** $< 5\%$
- Cd free atoms (dominant species) + molecular/aerosol-associated
- VSG-AAS LOD $60 \text{ pg mL}^{-1} \text{ Cd}$
- VSG-AFS LOD $0.42 \text{ pg mL}^{-1} \text{ Cd}$


VSG-AFS of Cd – METHOD VALIDATION

CRM	certified value (ng mL ⁻¹)	found value (ng mL ⁻¹)
NIST 1643 f	5.85 ± 0.13	5.90 ± 0.44
ERM-CA 713	5.00 ± 0.05	5.09 ± 0.20


- VSG of non-hydride forming element – Cd
- **Novel hydride atomizers for AAS based on DBD plasma**
- VSG-based speciation analysis
- VSG of mercury based on PMVG of Hg

HYDRIDE ATOMIZATION - AAS

- **externally heated quartz tube atomizers (QTA)**
 - heated to 900 °C
 - advanced construction (MMQTA) supplied by air/O₂
 - most common hydride atomizer



- **dielectric barrier discharge (DBD) plasmas**
 - low temperature, ambient pressure plasmas
 - AC high voltage
 - novel atomizer



ANALYTICAL FIGURES OF MERIT - SENSITIVITY

Atomizer	Sensitivity, s ng ⁻¹					
	Pb	Bi	Sn	Se	Te	As
(MM)QTA	0.29 ± 0.01	0.40 ± 0.02	0.33 ± 0.01	0.53 ± 0.03	0.32 ± 0.01	0.48 ± 0.01
DBD	0.09 ± 0.01	0.15 ± 0.01	0.05 ± 0.01	0.60 ± 0.04	0.32 ± 0.01	0.54 ± 0.04

- **As, Se, Te – sensitivity reached in DBD comparable to (MM)QTA**
- **Pb, Bi, Sn – (MM)QTA performs much better (3-7 times higher sensitivity)**
- **(MM)QTA – sensitivity difference among elements: factor of 2**
- **DBD - sensitivity difference among elements: factor of 12**


B. Baranová et al., Spectrochim. Acta B, accepted.

L. Juhászová et al., Spectrochim. Acta B 158 (2019), 105630.



P. Novák et al., Anal. Chem. 88 (2016), 6064-6070.

MECHANISTIC STUDIES - LIF

a) **Pb: atomization efficiency $23 \pm 7\%$**



b) **Te: atomization efficiency $100 \pm 7\%$**



Mechanistic studies – DEPOSITED FRACTION

leaching experiments, ICP-MS detection

Analyte	Analyte fraction (%) deposited in the atomizer	
	DBD	(MM)QTA
Pb	91 ± 5	107 ± 4
Bi	94 ± 1	92 ± 3
Se	26 - 43	15 ± 2
Te	62 ± 2	37 ± 2


- fast decay of Pb and Bi free atoms → deposit formation → low sensitivity (DBD)
- spatial distribution of deposits differ between DBD and (MM)QTA
 - DBD – homogeneous distribution even in the discharge area
 - MMQTA – in the colder atomizer zones

J. Kratzer et al., Anal. Chim. Acta 1028 (2018) 11-21.


J. Kratzer et al., Anal. Chem. 88 (2016) 1804-1811.

B. Baranová et al., Spectrochim. Acta B, accepted.

Spatial distribution of ^{75}Se in DBD and MMQTA after atomization of 3 replicates of a ^{75}Se tracer sample solution (exposure time of 68 hours).



PRECONCENTRATION *IN-SITU* IN DBD



1) ANALYTE TRAPPING


Ar + O_2

2) ANALYTE RELEASE

Ar + H_2 (blank)

No change in DBD HV / power settings

IN-SITU PRECONCENTRATION IN DBD



IN-SITU PRECONCENTRATION IN DBD

element	Preconcentration efficiency, %	LOD, ng mL ⁻¹
As	100	0.01
Se	70	0.01
Sb	100	0.02
Te	51	-
Bi	60	-

P. Novák et al., Anal. Chem. 88 (2016), 6064-6070.

J. Kratzer et al., J. Anal. Atom. Spectrom. 34 (2019), 193-202

P. Zurynková et al., Anal. Chim. Acta 1010 (2018) 11-29.

K. Bufková et al., Spectrochim. Acta B 171 (2020) 105947.

J. Kratzer et al., Anal. Chem. 86 (2014), 9620-9625.

- VSG of non-hydride forming element – Cd
- Novel hydride atomizers for AAS based on DBD plasma
- **VSG-based speciation analysis (Hg, Te, Ge, As)**
- VSG of mercury based on PMVG of Hg

A) selective VSG

VSG → detection

Te(IV) and Te(VI)

B) post-column VSG

separation → VSG → detection

HPLC-VSG-ICP/MS

C) Generation of substituted volatile species


VSG → separation → detection

VSG-CT-ICP/MS

selective VSG

VSG → *detection*

Te(IV) and Te(VI)



ICP-MS/MS

LOD: 2.3 µg L⁻¹

VSG-ICP-MS/MS

LOD: 0.07 ng L⁻¹

A) selective VSG

VSG → detection

Te(IV) and Te(VI)

B) post-column VSG

separation → VSG → detection

HPLC-VSG-ICP/MS

C) Generation of substituted volatile species

VSG → separation → detection

VSG-CT-ICP/MS

A) selective VSG

VSG → detection

Te(IV) and Te(VI)

B) post-column VSG

separation → VSG → detection

HPLC-VSG-ICP/MS


C) Generation of substituted volatile species

VSG → separation → detection

VSG-CT-ICP/MS

CRYOGENIC SEPARATION

TRAP




trapping

- Cryogenic trap (CT)
- U-tube immersed in liquid N₂
- all species retained @ -196 °C
- after heating the U-tube
- species separated according to their b.p.



RELEASE




release/separation




Speciation analysis of Ge – method development

- VSG → separation → detection



Speciation analysis of Ge – applications


- VSG → separation → detection



Speciation analysis of As – applications

- VSG → separation → detection

- in whole blood/plasma without extraction




- 50-100 μl samples
- LOD $\sim \text{pg ml}^{-1}$
- normal levels of exposure

Generation of alkyl-/aryl-substituted volatile species

- **VSG** → separation → detection
- **VSG** from HCl and TRIS buffer media
- **cryogenic trap (CT)** used for separation

- $\text{Hg}^{2+} \rightarrow \text{Hg}^0$
- $\text{MeHg}^+ \rightarrow \text{MeHgH}$
- $\text{EtHg}^+ \rightarrow \text{EtHgH}$
- $\text{PhHg}^+ \rightarrow \text{PhHgH}$



- **decomposition of substituted species during VSG step !!!**
- more pronounced in **HCl** than **TRIS** buffer media

quantification of fraction decomposed to Hg⁰ (%)

	HCl	TRIS
MeHgH	41	6
EtHgH	77	28
PhHgH	94	99

A) selective VSG

VSG → detection

Te(IV) and Te(VI)

B) post-column VSG

separation → VSG → detection

HPLC-VSG-ICP/MS



C) Generation of substituted volatile species

VSG → separation → detection

VSG-CT-ICP/MS

Post-column VSG

- separation → VSG → detection



HPLC-VSG-ICP/MS

- mobile phase matrix elimination
- constant IS response
- increased ICP/MS sensitivity

Fig. 3. Chromatograms of mercury species (^{202}Hg) mixed standard solution of Hg^{2+} , MeHg^+ , EtHg^+ and PhHg^+ containing $1 \mu\text{g L}^{-1}$ (as Hg) of each species, obtained without (A) and with postcolumn VSG step (B). The IS signal (^{125}Te signal) obtained without (C) and with postcolumn VSG step (D).

Post-column VSG:

separation → VSG → detection

Analytical figures of merit found for HPLC-ICP/MS and HPLC-VSG-ICP/MS


Species	Slope Cts L μg^{-1}	LOD (ng L $^{-1}$)
HPLC-ICP-MS		
Hg $^{2+}$	10 213	15
MeHg $^{+}$	9 957	15
EtHg $^{+}$	10 120	17
PhHg $^{+}$	3 499	26
HPLC-VSG-ICP-MS		
Hg $^{2+}$	398 430	3
MeHg $^{+}$	351 989	2
EtHg $^{+}$	336 402	4
PhHg $^{+}$	178 399	6

HPLC-VSG-ICP/MS

Sensitivity increased 30-40 times

LOD improved 3-7 times

VSG-based speciation analysis of Hg




Generation of substituted VS

analytical artifacts

M. Migašová et al., Anal. Chim. Acta 1119 (2020), 68-76.


post-column VSG


reliable approach



I. Petry-Podgórska et al., Microchem. J. 170 (2021), 106606.

- VSG of non-hydride forming element – Cd
- Novel hydride atomizers for AAS based on DBD plasma
- VSG-based speciation analysis
- **VSG of mercury based on PMVG of Hg**






laboratory made DBD reactor

AAS detector

QTA atomizer @ ambient temperature
 150 mL min⁻¹ He

Samples – droplets (2 μL)

Hg²⁺
 MeHg⁺



38 kV, 23 W, 40 kHz, 60% duty cycle

laboratory made power supply source

CONCLUSIONS

- **VSG of Cd**

- promising approach, 60% efficiency

- **novel DBD hydride atomizers**

- can compete with QTAs (As, Se, Te)
- *in-situ* preconcentration feasible

- **VSG for speciation analysis**

- postcolumn VSG – reliable approach
- generation of substituted VS – **artifacts** due to species decomposition (**Hg**)
 - **reliable** approach for **As** and **Ge**

- **PMVG of Hg**

- high introduction efficiency
- good choice for volume-limited samples

ACKNOWLEDGMENTS

Milan Svoboda



Czech Academy
of Sciences

Jiří Dědina

Gilberto da Silva Coelho Junior

Tomáš Matoušek

Inga Petry-Podgórska



FACULTY OF SCIENCE
Charles University

Barbora Baranová (Kodříková)

Zuzana Králová

Linda Sagapova

Věra Schrenková

Michaela Migašová



Joachim Franzke

Sebastian Burhenn

Sebastian Brandt



Pavel Dvořák

Martina Mrkvičková

Jan Voráč

Institute of Analytical Chemistry of the CAS

Department of Trace Element Analysis

ACKNOWLEDGMENTS

