


**M U N I**  
**M E D**

**Indirect calorimetry.**

Laboratory conditions:  
Room temperature (t) = 22.5 °C  
Barometric pressure (B) = 101.5 kPa



# Resting condition




Chart View - Krogh Spirometer x

No Sampling

A - Krogh Spirometer Averag... x

0.004514470756 L/s

Channel: 1 1 Comment Add 100 /s

# Standing position




Chart View - Krogh Spirometer  
No Sampling

A - Krogh Spirometer Averag...  
0.006449207771 L/s

No Sampling  
Krogh Spirometer

# After workload




Chart View - Krogh Spirometer x  
No Sampling  
A - Krogh Spirometer Averag... x  
0.01155942654 L/s