

day	group_nr	treatment	rat_mass_ξ	APTT	thrombus_mass
Po_ráno	1	heparin	268	200	0.015
	2	heparin	267	200	0.01
	4	heparin	282	200	0.015
	3	kontrola	255	17.4	0.03
	4	kontrola	261	16.5	0.04
	5	kontrola	279	99	0.025
po_dopo	1	heparin	315	200	0.01
	2	heparin	292	200	0.005
	3	heparin	278	200	0.01
	4	heparin	288	200	0.015
	4	kontrola	275	16.7	0.04
	5	kontrola	295	20.4	0.055
	3	kontrola	274	29.3	0.05
út_ráno	1	heparin	315	200	0.015
	3	kontrola	313	34.9	0.025
	4	kontrola	311	20.7	0.025
	5	kontrola	304	16	0.025
	3	heparin	269	200	0.01
	4	heparin	295	200	0.015
	2	kontrola	283	28.8	0.045
út_dopo	3	heparin	295	200	0.01
	4	heparin	300	200	0
	1	kontrola	300	43.2	0.025
	2	kontrola	300	17.6	0.025
	3	kontrola	317	16.3	0.035
	4	kontrola	326	21.3	0.03
	5	kontrola	302	96.1	0.025
čt_rano	1	heparin	315	200	0.01
	2	heparin	302	200	0.005
	3	heparin	247	200	0.01
	4	heparin	274	200	0.01
	1	kontrola	287	21.8	0.05
	2	kontrola	285	50.4	0.025
	3	kontrola	277	20.6	0.08
čt_dopo	1	heparin	295	200	0.01
	3	heparin	300	200	0.015
	4	heparin	278	200	0.01
	1	kontrola	280	17.1	0.025
	2	kontrola		46.2	0.025
	3	kontrola	285	14.7	0.05
	5	kontrola	290	13.7	0.06
pá_ráno	2	heparin	309	200	0.015
	1	kontrola	305	19.7	0.04
	2	kontrola	282	37.4	0.04
	3	kontrola	305	46.2	0.04

	4 kontrola	304	13.5	0.05
	5 kontrola	289	36.1	0.02
pá_dopo	1 kontrola	259	19.3	0.035
	2 kontrola	265	33.6	0.015
	3 heparin	258	200	0.015
	4 kontrola	277	14	0.03
	5 kontrola	255	11.2	0.045
	2 heparin	258	200	0.005
	1 heparin	254	200	0.01
	5 kontrola	258	11.2	0.045
po_ráno	5 kontrola	260	18.5	0.045
	2 heparin	278	200	0.015
	3 heparin	305	200	0.01
	3 kontrola	292	77.4	0.05
	2 kontrola	283	36.9	0.045
	4 kontrola	295	20.8	0.035
	1 heparin	300	200	0.01
	1 kontrola	285	45.7	0.045
	4 heparin	300	200	0.01
po_dopo	4 kontrola	346	25.3	0.035
	2 kontrola	281	35.9	0.02
	3 kontrola	300	34.5	0.02
	5 kontrola	325	27.3	0.025
	1 heparin	312	200	0.01
	2 heparin	305	200	0.01
	3 heparin	282	200	0.015
UT samoplak	1 heparin	300	200	0.005
	2 heparin	300	200	0.01
	3 heparin	300	200	0.015
	4 heparin	300	200	0.01
	1 kontrola	300	58.3	0.04
	2 kontrola	300	25.2	0.045
	3 kontrola	300	11.3	0.03
ST rano	1 kontrola	290	45.5	0.055
	1 heparin	305	200	0.015
	2 kontrola	325	23.3	0.025
	4 heparin	282	200	0.015
	5 kontrola	270	21.2	0.03
	3 heparin	310	200	0.01
	3 kontrola	310	15.2	0.035
	4 kontrola	295	15.2	0.035
Wednesday	1 kontrola		13.5	0.02
	2 heparin	315	200	0.015
	2 kontrola		26.5	0.025
	3 heparin		200	0.01
	3 kontrola	315	200	0.025
	4 heparin	289	200	0.01
Pá_ráno	1 kontrola	325	25.3	0.045
	2 kontrola	297	18.9	0.05
	3 kontrola	281	29.4	0.045

	4 kontrola	288	33.2	0.05
	5 kontrola	277	19.9	0.04
	1 heparin	317	200	0.01
	2 heparin	279	200	0.01
	4 heparin	283	200	0.01
Pa_dopo	4 heparin	280		0.01
	3 heparin	317		0.01
	1 kontrola	287	15.2	0.055
	5 kontrola	265	33	0.045
	3 kontrola	290	22	0.035
	4 kontrola	287	25.4	0.05