


HRANIČNÍ PLAZMOLÝZA

c [M]	česnek cibule		douška hustoli	
	měř. 1	2	měř. 1	2
0	0	0	0	0
0.1	0	0	0	0
0.2	0	0	0	0
0.3	0	0	65	70
0.4	0	0	100	100
0.5	0	0	100	100
0.6	90	100	100	100
0.7	100	100	100	100
0.8	100	100	100	100
0.9	100	100	100	100

REFRAKTOMETRIE n vzorek (lilek brambor)

c [M]	n vzorek (lilek brambor)		
	n kal.	1	2
0	1.9	1.8	1.8
0.1	4.9	5.2	5.2
0.2	8	8.3	8.2
0.3	12	11.6	11.6
0.4	14.5	13.4	13.6
0.5	17.8	16	16
0.6	20.9	17.3	17.2
0.7	24	19.8	20
0.8	27.2	21.1	21.4
0.9	30	23.8	23.8



TLAKOVÁ METODA

t [min]	tlak [MPa]			zalitý	1	2	3
	suchý	1	2				
0	-0.36	-0.3	-0.295	-0.25	-0.41	-0.47	
15	-0.46	-0.31	-0.32	-0.31	-0.48	-0.56	

stá

t = 24 °C

