

Spektrofotometrické stanovení obsahu fotosyntetických pigmentů

	<b>663 nm</b>	<b>646 nm</b>	<b>470 nm</b>	<b>V (ml)</b>	<b>m (g)</b>
Kontrola	1.452	0.548	0.973	25	0.25
Kontrola	1.803	0.676	1.245	30	0.25
Kontrola	0.691	0.28	0.519	25	0.25
bez P	1.065	0.453	0.776	25	0.25
bez P	1.193	0.552	1.021	25	0.27
bez P	1.181	0.46	0.829	25	0.23
bez N	0.622	0.253	0.521	25	0.25
bez N	0.725	0.269	0.553	30	0.27
bez N	0.448	0.245	0.415	25	0.24
bez Fe	0.183	0.067	0.242	25	0.23
bez Fe	0.21	0.082	0.309	25	0.27
bez Fe	0.143	0.059	0.151	25	0.24

**LA (cm2)**

- 8.59572
- 19.73832
- 20.33
- 11.5974
- 19.95056
- 17
- 12.886
- 19.278972
- 18.23
- 10.4604
- 17.045525
- 14.74