


cellSens Entry

File Edit View Acquire Image Process Measure Tools Window Help

56,3 % 4x 10x 20x 40x Reflected Live Snap 300,8 ms

Start Page Live (active)



50 100 150 200 250 300 350 400 450 500 550 600 650 700 μm

50 100 150 200 250 300 350 400 450 500 μm

50 μm

Measurement

Camera Control

Live Snap

Movie recording

Exposure: 300,8 ms

Manual Automatic

300,8 ms

Exposure time:

Gain: 0 dB

Resolution: 2080 x 1544

Live/Movie: 2080 x 1544

Snap: 2080 x 1544

Apply to Live

Note

Camera Con... Adjust Display

Ready 56,3 % (5,63x) 9:10 14.4.2023