



	Skupina 1		Skupina 2		Skupina 3		Skupina 4	
M	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I



	Skupina 5		Skupina 6		Skupina 7		Skupina 8	
M	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I	PCR	<i>Hinf</i> I




Skupina 9

M

PCR

*Hinf*I



Délky Markeru:

100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 2000, 3000 bp